


visicolor[®] ECO Nitrat


Nitrate / Nitrate / Nitrato

Method(e) / Método	PF-3 Soil 450 nm	PF-3 Fish 450 nm	PF-12 436 nm	PF-12 ^{Plus} 436 nm	Advance 450 nm
5411 1.0–14.0 mg/L NO ₃ -N	■	■	■	■	■
5412 4.0–60.0 mg/L NO ₃ ⁻	■	■	■	■	■




Null mit Probe
Blank with sample
Blanc avec échantillon
En blanco con muestra


Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar




Küvetten spülen
Rinse test tube
Rincer les éprouvettes
Lavar los tubos




Probe
Sample
Echantillon
Muestra




5 δ NO₃-1



Schütteln
Shake
Agiter
Agitar



1 — NO₃-2




Schütteln
Shake
Agiter
Agitar



Säubern
Clean
Nettoyer
Limpiar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición