

## Data sheet

## Filter paper

| Parameter*                    | Value*                                          |
|-------------------------------|-------------------------------------------------|
| Grade                         | MN 714                                          |
| Basis weight                  | 75 g/m <sup>2</sup>                             |
| Thickness                     | 0.25 mm / 9.8 mil                               |
| Filtration speed              | 20 s                                            |
| Surface                       | Embossed                                        |
| Color                         | White                                           |
| Material                      | Cellulose                                       |
| Characteristics / Application | Medium fast filtration, e.g. for essential oils |

\*These data are typical values and must not be regarded as specification.