

## Data sheet

## Filter paper

Parameter*	Value*
Grade	MN 616 md
Basis weight	85 g/m <sup>2</sup>
Thickness	0.18 mm / 7.1 mil
Filtration speed	55 s
Surface	Smooth
Color	White
Material	Cellulose
Ash content	ca. 0.1 %
Average retention capacity	4–7 µm
Characteristics / Application	Medium fast filtration, e.g. for magnesium-ammonium-phosphate, zinc-sulphide etc.

\*These data are typical values and must not be regarded as specification.