

**NANOCOLOR®** gesamt-Kjeldahl-Stickstoff TKN 16  
 total Kjeldahl nitrogen TKN / Azote total Kjeldahl TKN / Nitrógeno total Kjeldahl NTK



Test 0-67  
 REF 985067

Methode / Method / Méthode / Método 345 nm  
 0671 1.00–32.0 mg/L TN<sub>b</sub>  
 0672 1.00–16.0 mg/L NO<sub>3</sub>-N  
 0673 1.00–16.0 mg/L TKN

Min	20–80 %		Max
1.00	7.20	25.8	32.0 mg/L TN <sub>b</sub>
1.00	4.00	13.0	16.0 mg/L NO <sub>3</sub> -N
1.00	4.00	13.0	16.0 mg/L TKN

