

MACHEREY-NAGEL

Urine test strips Combi 10[®] SGL



Urine analysis

- Precise
- Cost-effective
- No interference from vitamin C

Medi-Test

MACHEREY-NAGEL


www.mn-net.com



Urine test strips Combi 10[®] SGL

Best available vitamin C protection

Comparative data clearly show: Combi 10[®] SGL urine test strips have the best protection of blood and glucose tests against vitamin C. With other brands, vitamin C may prevent the oxidation reaction essential for these test. This may cause false negative results for blood and glucose. Using Medi-Test Combi 10[®] SGL, blood and glucose in urine will be safely detected. Patients can receive treatment at the earliest stage thus increasing chances for success.



Concentration of vitamin C in the sample with no effect on the results for blood and glucose in urine (mg/dL)


Desiccant enclosed in cap

Test strips Medi-Test Combi 10[®] SGL come in tubes with a desiccant containing cap. The desiccant ensures a low humidity in the tube and the longevity of the test strips. Many other packages contain a simple desiccant pack that may easily be spoiled or accidentally removed. The desiccant in the cap is convenient and safe. Closing the cap safely protects the test strips.



High quality color chart

The reaction colors on the test strip exactly match the colors on the label. This is ensured by a high quality printing process using specially mixed full colors. Other color charts often yield far less precise matches that may lead to inaccurate results. The evaluation of Medi-Test Combi 10[®] SGL is always clear and reliable providing perfectly accurate results.



Presented by

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

