



visicolor® ECO Eisen 1

Iron / Fer / Hierro

Messbereich / Range / Domaine de mesure / Rango
0.04 – 1.0 mg/L Fe



Probe
Sample
Echantillon
Muestra




5 Δ Fe-1



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar



MACHEREY-NAGEL
www.mn-net.com