



visocolor[®] ECO pH 4.0–9.0
pH 4.0–9.0 / pH 4.0–9.0 / pH 4.0–9.0

Messbereich / Range / Domaine de mesure / Rango
pH 4.0–9.0




Probe
Sample
Echantillon
Muestra



4 Ø pH-1



Schütteln
Shake
Agiter
Agitar



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar

