



visocolor® ECO Hydrazin

Hydrazine / Hydrazine / Idrazin


Messbereich / Range / Domaine de mesure / Rango
0.05 – 0.40 mg/L N₂H₄




filtrieren wenn trüb
filter if cloudy
filtrer si trouble
filtra si turbido




Probe
Sample
Echantillon
Muestra




5 ∅ N₂H₄-1



Schütteln
Shake
Agiter
Agitar




5 ∅ N₂H₄-2



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar



MACHEREY-NAGEL
www.mn-net.com

09.2022