

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Nickel

| Parameter* | Value* |
|------------------------------|---|
| REF | 91305 |
| Evaluable on reflectometer | No |
| Shelf life (from production) | 2.5 Years |
| Scope of delivery | 100 test strips in a tube, instruction leaflet |
| Parameter | Nickel |
| Gradation | 0 • 10 • 25 • 50 • 100 • 250 • 500 • 1000 mg/L Ni ²⁺ |
| Method | Dimethylglyoxime |
| Color change | white → bright-red |

*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11
52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com