

MACHEREY-NAGEL

Vials and caps



For reliable analysis



MACHEREY-NAGEL

www.mn-net.com



Quality since 1911

Since 1911 MACHEREY-NAGEL and their products stand for high quality, innovation and reliability in chemical and bio-molecular analysis. Nowadays we offer an extensive product portfolio in the areas of filtration, rapid tests, water analysis, chromatography and bio analysis. MACHEREY-NAGEL is a family-owned company having the mission to fulfill the versatile requirements of our customers in an optimal way. Personal and competent advice as well as customer-friendly support belong to our philosophy as does the high quality of our products.

With approx. 700 employees and 40 labs we serve our customers all around the world, either directly or indirectly through our distribution channels in more than 150 countries (please see www.mn-net.com/distributors). Our main head-quarter is located in Düren (Germany) with subsidiaries in France, Switzerland and in the USA.




Filtration



Rapid Tests



Water Analysis



Chromatography



Bioanalysis

Chromatography with MACHEREY-NAGEL – Complete solutions for analysis

MACHEREY-NAGEL has a long-term expertise in chromatography, as we have been involved in this business right from its beginnings. Today we are a supplier for all kinds of chromatography consumables, from sample preparation to final analysis. Our product range includes HPLC and GC columns as well as syringe filters, SPE cartridges or TLC plates.

Of course we also offer a comprehensive delivery program of autosampler vials and closures along with some accessories, like e.g. crimping tools.



HPLC



GC



TLC



SPE and Flash



Syringe filters



Vials and caps

Vials and caps

MACHEREY-NAGEL offers a wide product range of vials, caps and accessories for chromatography. You can find for almost all autosamplers and applications the right product in our delivery program. In order to facilitate the correct selection of products, we have implemented some useful tools in the catalog (pls. see Autosampler Compatibility Charts starting on page 8 and Real Size Drawings starting on page 57) as well as on our website (pls. see www.mn-net.com/vialfinder).

All vials are made from 1st hydrolytic class glass except if explicitly stated otherwise. All septa material used for the production of closures or septa has – aside from physical inspection – also been tested analytically. Certificates of conformity for all glass products and closures can be downloaded from our website under www.mn-net.com/certificate or be requested at support@mn-net.com Cost-free samples of all vials and caps can be asked for at any time.

MACHEREY-NAGEL – About us

MACHEREY-NAGEL Vials and caps for



GC and HPLC



Headspace

Production

All our products are manufactured under strict hygienic conditions. Opto-electronical supervision at the glass processing machines ensures high quality for a safe analysis. Of course additional in-process controls and a final QC inspection are also part of the measurements which we have implemented to grant a reliable product quality.

Fully automated closure manufacturing with defined production parameters avoids variations in the physical and analytical condition of the most critical component of the analysis: the cap. Strict analytical specifications on the septa material allow only processing of such septa material that passes the comprehensive analytical tests in our lab.

Almost all of our closures, caps and plastic vials are automatically counted (not weighed) for exact quantity obedience and tamper-evidently packed in zip-lock bags. The appropriate product label is directly printed onto the bag's foil.



Environmental analysis (TOC/VOC)



Sample storage



Fully automated production for safe and reliable consumables for chromatography



Packaging

Vials and caps can be bought as individual components with 100 pieces in a box respectively in a PE bag. All packaging material we use is also strictly evaluated with regard to possible releases of contaminants to ensure that the products remain pure during storage for a safe and reliable analysis. We also took care that the zip lock bags for the closures as well as the boxes for the vials are resealable, so that consumers can protect the products against any contamination in the lab during consumption.


Pre-sealed vial-closure combinations

Pre-sealed vial-closure combinations of any type are available as well. You can choose any type of vial with a closure of your choice and ask for your individual quotation for a pre-sealed combination of these.

Pre-sealing is done in-house under strict hygienic conditions.

Certificates

In case you need a certificate of conformity (CoC) for your ISO documentation, you can either download it from our website under www.mn-net.com/certificate or contact support@mn-net.com for their assistance. CoCs are based on our own internal QC reports and are batch related. In order to give you an impression on how they look like, please see below an example of such a CoC.



MACHEREY-NAGEL

Certificate of Analysis

Product information

Product	1.5 mL Screw Neck Vial N 9
	outer diameter: 11.6 mm
	outer height: 32 mm
	amber, flat bottom, wide opening, label + scale
REF	702284
LOT	2023017876
Date of examination:	23.01.2024

Confirmation

MACHEREY-NAGEL certifies that the above mentioned products have been manufactured against standards according to the valid specification, drawing and agreed lists of defects and that the specification limits have been met. The components have been manufactured with fully approved tooling. All required inspections and tests have been performed and the items delivered conform to the specifications of the order.

MACHEREY-NAGEL furthermore confirms that the vials / inserts have been manufactured from 1st hydrolytic class glass which meets the requirements of analytical purity and chemical resistance needed when used in chromatography applications.

MACHEREY-NAGEL is certified to EN ISO 9001:2015. A copy of this certificate can be provided upon request.

This document has been produced electronically and is valid without a signature.

MACHEREY-NAGEL GmbH & Co. KG
Valenciennestraße 11
52355 Düren - Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 66 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com

Contents

New products	6
Autosampler compatibility	8
 Shell vials N 8 / Crimp neck vials and caps N 8	18
 Screw neck vials, inserts and caps N 8	19
 Screw neck vials, inserts and caps N 9	21
 Screw neck vials, inserts and caps N 10	25
 Crimp neck vials, inserts and caps N 11	26
 Snap ring vials, inserts and caps N 11	30
 Crimp neck vials and caps N 13	32
 Screw neck vials, insert and caps N 13	33
 Diverse vials and caps	34
 Screw neck vials and closed top screw closures for storage of liquid samples	34
 Snap cap vials (for storage of powdery samples) N 18 + N 22 with PE snap caps ...	36
 Headspace screw neck vials and magnetic screw caps N 18 ...	37
 Crimp neck vials and caps N 20	38
 Screw neck vials and caps N 24	43
 Vial accessories and other consumables for chromatography ...	46
 GC injection port septa	46
 Capillary connectors for GC	46
 Manual crimping tools - Standard	47
 Manual crimping tools - Ergonomic	48
 Electronic crimping tools - Battery-powered	48
 Electronic crimping tools - High Power crimping tool	49
 Accessories for crimping tools	49
 Vial containers	50
 Vial racks	50
Numerical index	51
Practical tips for crimping	56
Real size drawings	57
General remarks	62

New products

Headspace screw neck vials and magnetic screw closures N 18 for CTC Combi PAL® and equivalent autosamplers

- New vial types, like 20 mL clear vial with label and scale or 10 mL amber vial
- Open and closed top screw closures N 18 with different septa materials are available
- Aside from Silicone/PTFE and Red Rubber/TEF, also closures with Butyl/PTFE are for sale
- Replacement septa of almost all closures can be obtained
- Magnetic screw closures N 9 and magnetic crimp closures N 11 for 1.5 mL vials as well as magnetic or bi-metal crimp closures N 20 are also included in our delivery program



Vials with firmly integrated inserts

Firmly integrated micro-inserts have the advantage of a ready-to-use unit that can be filled easily and safely due to the firm fixation of the insert and that can directly be placed on the lab bench or in a tray/rack afterwards. Sample removal is carried out without any risk for the needle due to the reliable centering of the micro-insert. The risk of any contamination is always lower with ready-to-use sample containers in comparison to individual components that still need to be assembled in the lab.

- Available as glass-in-glass product with insert firmly being melted to the outer glass vial in clear and amber. Clear vials with writing patch have been added to the delivery program.
- Available as glass in plastic product with a glass insert firmly pressed into an outer plastic vial shell for high robustness and excellent price-performance ratio. Upon request also available with integrated silanized glass inserts.



New natural rubber/ETFE closures N 11

- Available with silver (70246), green (70246.1), red (70246.2) and blue cap (70246.3)
- Very cost-effective closure for routine analysis in GC and HPLC
- Excellent resealability properties
- Temperature resistant from -40°C up to 120°C
- Alternative product for 702001, 70256, 70231.1, 70231.2, 70231.3



New products

Silicone/PTFE Pharma-Fix septa

Headspace analysis requires closures with a high temperature resistance, outstanding sealing properties and a high analytical purity. Our new Pharma-Fix Silicone/PTFE septa fulfill all these demands with a temperature resistance of - 60°C up to 200°C and the proven, outstanding performance of Silicone/PTFE with regard to cleanliness and penetrability. With these so-called Pharma-Fix septa the PTFE lamination is only applied centrally in the contact area with the sample, so that the uncoated, elastomeric septa rim can seal the vial's glass crimp neck tighter than completely PTFE laminated septa with their slippery surface. The side edges of these liners are designed lamellar so that they can be well compressed during the crimping process and be pressed firmly onto the glass rim. The drawn-in edge of the center hole keeps the liner in a stable, flat position even in case of high internal vial pressure.

Get the optimal closure for your headspace analysis that combines the following properties:

- Outstanding sealing properties due to only centrally applied PTFE (Pharma-Fix)
- High analytical purity due to long approved silicone/PTFE septa material
- Temperature resistance from - 60 °C up to 200 °C
- Safe and easy to penetrate due to the softness of Silicone/PTFE
- Lamellar designed side edges for good compressing during crimping
- The drawn-in edge of the center hole keeps the liner in a stable, flat position



PTFE-free (fluorine-free) closures for PFAS analysis


If you are going to carry out a PFAS analysis, it is decisive to select the right vials and closures for this application. Adsorption effects of glass as well as a possible contamination of the sample by particles from the septa, especially from the PTFE lamination, can endanger your analysis. Polyimide or polypropylene coated, fluorine-free septa ensure reliable analysis results and good chemical resistance. Available as 9 mm screw caps (MN REF 702402, 702405 and 702406), 11 mm snap ring cap (MN REF 702403) and as 20 mm bimetal crimp cap (MN REF 702404).



Snap ring closures N 11 with different cap colors

Colors play an important role in our lives and are perceived to mean different things. This can be valid, too, in the working life in a lab. Color-coded samples may represent different types of samples, samples from different shifts or from different labs and help us to differentiate them from each other. For snap ring caps N 11 we now offer a broad variety of cap colors with different septa materials to allow color-coding. Aside from this practical aspect, differently colored caps may also have an unconsciously influence on your emotions, as colors make life brighter and people happy.

In addition we also offer a light blue, soft cap with different liners inside that is easier to push on and to remove.



Autosampler compatibility

Agilent®

Main chapter	Pages	Most popular MN products for use on Agilent® instruments	
GC		5880, 5890, 6850, 6890, 7650ALS, 7673A, 7683A, 7693A, 7696A, CTC GC PAL® (G650X Agilent® GC Injector 80/120), PAL® 3 (LSI, RSI, RTC)	
 N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702878, 702025
 N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702006, 702088, 702007, 702007.0, 702008, 702135, 702335, 702009, 702172, 702010, 702500
	21	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819, 702825
	22, 23	Closures:	702732, 702080, 702082, 702081, 702287.1, 702037, 702038, 702402 (PFAS), 702406 (PFAS), 702035, 702084, 702085, 702026, 702155 (for GC PAL®)
 N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702075, 702892, 702076, 702888, 702015, 702016, 702891, 702891.0, 702014, 702134, 702334, 702809, 702173, 702174
	26	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819, 702825
	28	Closures:	70256, 702001, 70246, 70246.3, 70231.1, 70246.1, 70231.2, 70246.2, 702730, 702730.3, 702730.1, 702730.2, 70288, 70288.1, 70288.2, 70288.3, 702995, 702995.1, 702995.2, 702995.3, 702879 (for GC PAL®)
 N 13 Screw (only 6850, 22 Pos. Tray + 7693A)	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702050
HPLC		1050 Series (79855A), 1090 Series, 1100 Series (G1313A/G1329A/G1389A/G2260A), 1200 Series (G1367B, G1367D, G1377A), 1260 Infinity Series, 1290 Infinity Series, 1290 Infinity II Multisampler (G7161B), CTC HTS/HTC PAL®	
 N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702878, 702025
 N 9 Screw (standard samples)	21	Vials:	As indicated under GC – Screw N 9 – Vials
	21	Inserts:	As indicated under GC – Screw N 9 – Inserts
	22, 23	Closures:	702732, 702080, 702082, 702081, 702287.1, 702037, 702038, 702402 (PFAS), 702405 (PFAS), 702406 (PFAS), 702084, 702085, 702035, 702040, 702083, 702288.1, 702026, 702027
 N 11 Crimp (standard samples)	26, 27	Vials:	As indicated under GC – Crimp N 11 – Vials
	26	Inserts:	As indicated under GC – Screw N 9 – Inserts
	28	Closures:	As indicated under GC – Crimp N 11 – Closures, however, not 702879
 N 11 Snap (standard samples)	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
	30	Inserts:	As indicated under GC – Screw N 9 – Inserts
	31	Closures:	702731, 702063, 702063.2080, 702295, 702710.1, 702710.2080, 702095, 702142, 702108, 702064, 702717.2, 702717.2080, 702143, 702150, 702151, 702403 (PFAS), 702718, 702718.1, 702718.2080, 702401
Headspace		HS7694, 8697, HS7697A, G1888A, CTC Combi PAL®, PAL® 3 (LSI, RSI, RTC)	
 N 18 Screw (Combi PAL®, G1888A, HS7697A, 8697)	37	Vials:	702866, 702866.2, 702826, 702826.0, 702826.2
	37	Closures:	702055
 N 20 Crimp	38	Vials:	702918, 702261 702924, 702263, 702263.0, 702263.2 (all for Combi PAL®, PAL® 3 LSI, RSI, RTC, HS7697, 8697)
	40, 41, 42	Closures:	70234.10, 70234, 702071, 702094, 702340, 702835, 702931.3 + 702804 = 70237.2 702929, 702341 (both for Combi PAL®, PAL® 3 LSI, RSI, RTC)
Purge + Trap (TOC)		AQUATEk 70, AQUATEk 100, SOLATEk 72, Archon Purge + Trap (7695A), ATOMX XYZ, PAL® 3 (LSI, RSI, RTC)	
 N 24 Screw	43	Vials:	702023, 702024
	44	Closures:	702058, 702060, 702062, 702171

Autosampler compatibility

CTC

Main chapter	Pages	Most popular MN products for use on CTC instruments	
GC		CTC (Leap) GC PAL®, CTC (Leap) Combi PAL®/ Combi-xt® PAL, CTC (Leap) PAL® 3 (RTC/RSI/LSI) as well as the appropriate OEM autosamplers (Gerstel, Shimadzu®, etc.) and equivalent instruments of other manufacturers (Konik)	
 N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702878
 N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702006, 702007, 702007.0, 702008, 702135, 702335, 702088, 702500, 702009, 702172, 702010
	21	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702825, 702819
	22, 23	Closures:	702155, 702156 (closures for liquid injection → please see HPLC)
 N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702892, 702075, 702076, 702888, 702891, 702891.0, 702014, 702134, 702134.1, 702334, 702015, 702016, 702809, 702173, 702174
	26	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702825, 702819
	28	Closures:	702879 (closures for liquid injection → please see HPLC)
HPLC		CTC (Leap) LC PAL®, HTX/HTX-xt PAL®, HTC /HTC-xt PAL®, HTS/HTS-xt PAL®, PAL® 3 (LSI, RSI, RTC), Dual Head PAL® as well as equivalent instruments of other manufacturers	
 N 9 Screw (standard samples)	21	Vials:	As indicated under GC – Screw N 9 – Vials
	21	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702825, 702819
	22, 23	Closures:	702287, 702287.1, 702036, 702037, 702038, 702107, 702026, 702402 (PFAS), 702405 (PFAS), 702406 (PFAS), 702288, 702288.1, 702039, 702040, 702083, 702109, 702027, 702286, 702035, 702158, 702084, 702085, 702159
 N 11 Crimp (standard samples)	26, 27	Vials:	As indicated under GC – Crimp N 11 – Vials
	26	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702825, 702819
	28	Closures:	70288, 70288.1, 70288.2, 70288.3, 702995, 702995.1, 702995.2, 702995.3, 702823
 N 11 Snap (standard samples)	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
	30	Inserts:	As indicated under GC – Screw N 9 – Inserts
	31	Closures:	702710.1, 702717.2, 702718.1, 702710, 702064, 702718, 702710.2080, 702403 (PFAS), 702717.2080, 702718.2080, 702095, 702142, 702108, 702143, 702150, 702151
Headspace		CTC (Leap) Combi PAL® / Combi-xt PAL®, PAL® 3 (LSI, RSI, RTC) as well as the appropriate OEM autosamplers (Gerstel, Shimadzu, etc.) and equivalent instruments of other manufacturers (Konik)	
 N 18 Screw	37	Vials:	702866, 702866.2, 702826, 702826.0, 702826.2
	37	Closures:	702827, 702055
 N 20 Crimp	38	Vials:	702924, 702263, 702263.0, 702263.2
	40, 41, 42	Closures:	702929, 702834, 702341
			For washer bottle 702924: 70267 + 702144
Purge + Trap (TOC)		PAL® 3 (LSI, RSI, RTC)	
 N 24 Screw	43	Vials:	702023, 702024
	44	Closures:	702058, 702060, 702062, 702171

Autosampler compatibility

Dionex® (Thermo Scientific®)

Main chapter	Pages	Most popular MN products for use on Dionex® instruments	
HPLC		UltiMate 3000 Basic System (ACC-3000, ACC-3000T), UltiMate 3000 Standard SD (WPS-3000TSL, WPS-3000SL), UltiMate 3000 Rapid Separation RS (WPS-3000RS, WPS-3000TRS), Biocompatible Systems (WPS-3000TBPL, WPS-3000TBRS), UltiMate 3000 RSLCnano (WPS-3000TPL RS, WPS-3000PL RS), UltiMate 3000 Semi-Preparative (WPS-3000TSL, WPS-3000SL), Autosampler with Fraction Collection (WPS-3000TFC, WPS-3000TBFC), Vanquish Systems (VC-A12, VC-A13, VF-A10, VF-A40, VH-A10, VH-A40), Summit ASI-100/ASI-100T, Autoselect AS50, Gina 50/160, Famos	
IC*		Dionex® AS, AS-AP, AS 40, AS 50, AS-HV	
■ N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	702025, 70289
■ N 8 Screw (standard samples)	19	Vials:	70213, 70213.2, 702004, 702893, 702860
	19	Inserts:	702968, 702968.1, 702824, 702005
	20	Closures:	70245, 702437
■ N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702007, 702007.0, 702008, 702135, 702335, 702006, 702088, 702009*, 702172*, 702010*, 702500*
	21	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819*, 702825
	22, 23	Closures:	702287.1, 702287, 702036, 702037, 702038, 702107, 702026, 702402 (PFAS), 702405 (PFAS), 702406 (PFAS), 702288.1*, 702288*, 702039*, 702040*, 702083*, 702109*, 702027*
■ N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702892, 702075, 702076, 702888, 702015, 702016, 702891, 702891.0, 702014, 702134, 702134.1, 702334
	26	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819*, 702825
	28	Closures:	70288, 70288.1, 70288.2, 70288.3, 702823, 70256, 70231.1, 70231.2, 70231.3
■ N 11 Snap (standard samples)	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
	30	Inserts:	As indicated under <i>Crimp N 11 – Inserts</i>
	31	Closures:	702710.1, 702710.2080, 702710, 702095, 702142, 702108, 702717.2, 702717.2080, 702064, 702143, 702150, 702151, 702403 (PFAS)
■ N 13 Screw (large sample volumes)	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702925
■ N 18 Screw (large sample volumes)	37	Vials:	702866, 702866.2
	37	Closures:	702827
■ N 20 Crimp (large sample volumes)	38	Vials:	702924 (not for Summit ASI-100/ASI-100T, Gina 60)
	40	Closures:	702094, 702093
■ N 24 Screw (large sample volumes)	43	Vials:	702023*, 702024*, 702074*, 702131* (for AS-HV)
	44	Closures:	702058*, 702171* (for AS-HV)
Purge + Trap (TOC, ASE)		ASE 200, ASE 350	
■ N 24 Screw	43	Vials:	702023, 702024 (only for ASE 200), 702074, 702131
	44	Closures:	702058, 702171

Autosampler compatibility

Knauer

Main chapter	Pages	Most popular MN products for use on Knauer instruments	
HPLC		TRIATHLON, Knauer Smartline 3800/3950, Knauer UHPLC Version PLATINblue AS 1, Azura® AS 6.1L, Knauer Optimas	
N 8 Screw (standard samples)	19	Vials:	70213, 70213.2, 702004, 702893
	19	Inserts:	702968, 702968.1, 702824, 702005
	20	Closures:	702067, 70245
N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702006, 702007, 702007.0, 702008, 702135, 702335, 702088, 702009, 702172, 702010, 702500
	21	Inserts:	702813, 702077, 702716, 702818, 702818.1, 702819, 702825
	22, 23	Closures:	702732, 702030, 702080, 702081, 702082, 702147, 702287.1, 702287, 702036, 702037, 702038, 702107, 702402 (PFAS), 702406 (PFAS), 702028, 702026
N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702892, 702075, 702076, 702888, 702891, 702891.0, 702014, 702134, 702134.1, 702334
	26	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819, 702825
	28	Closures:	70288, 70288.1, 70288.2, 70288.3, 702823, 702730, 702730.1, 702730.2, 702730.3, 70256, 70246.1, 70246.2, 70246.3
N 13 Screw (large sample volumes) (only for TRIATHLON)	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702050, 702925
N 18 Screw (large sample volumes)	37	Vials:	702866, 702866.2
	37	Closures:	702072, 702055, 702827
N 20 Crimp (large sample volumes)	38	Vials:	702918
	40, 42	Closures:	702094, 702129

PerkinElmer®

Main chapter	Pages	Most popular MN products for use on PerkinElmer® instruments	
GC		AS 100/100B, AS 300, AS 2000/2000B, AS8300, AI-1, AutoSytem, Turbomatrix MultiPrep / MultiPrep+, Clarus 400/480, Clarus 500/580/590, Clarus 600/680/690 (with 108-Vial Tray N6550017 for none PE-Vials, screw neck vials N 9 or N 10 may be used as well → see HPLC), AS2400 Liquid Sampler	
N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	702025, 70289, 702878
N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702078, 702079, 702284, 702006, 702009, 702172, 702010, 702500, 702007, 702007.0, 702008, 702135, 702335
	21	Inserts:	702818, 702818.1, 702819, 702825
	22, 23	Closures:	702732, 702028, 702030, 702080, 702081, 702082, 702147, 702287.1, 702026, 702287, 702036, 702037, 702038, 702107, 702402 (PFAS), 702406 (PFAS), 702155 (for Turbomatrix Multiprep/MultiPrep+)
N 11 Crimp (standard samples) * small opening ** wide opening *** not for Turbomatrix MultiPrep / Multiprep+	26	Vials:*	70201CG, 70214CG
	26	Inserts:*	702824, 702005
	26, 27	Vials:**	70201HP, 70201HP.2, 702885, 702892, 702075, 702076, 702888, 702891, 702891.0, 702014, 702134, 702134.1, 702334
	26	Inserts:**	702818, 702818.1, 702819, 702825
	28	Closures:	702730***, 70256***, 702730.1***, 702730.2***, 702730.3***, 70288, 70288.1, 70288.2, 70288.3, 702995, 702995.1, 702995.2, 702995.3, 702879 (for Turbomatrix Multiprep/MultiPrep+)
N 13 Screw	33	Vials:	702962, 702973
	33	Inserts:	702972, 702974
	33	Closures:	702926

Autosampler compatibility

PerkinElmer®

Main chapter	Pages		Most popular MN products for use on PerkinElmer® instruments
HPLC			Series 200 LC, ISS 100, ISS 200, ISS 225, LC 600, Integral 4000, Altus™ UPLC, Flexar™, LC 300
 N 8 Crimp (microsampling)	18	Vials:	70286
	18	Closures:	702025
 N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702078, 702079, 702284, 702009, 702172, 702010, 702500, 702007, 702007.0, 702008, 702135, 702335
	21	Inserts:	702818, 702818.1, 702819, 702825
	22, 23	Closures:	702288.1, 702027, 702288, 702039, 702040, 702083, 702109, 702287.1, 702026, 702287, 702036, 702037, 702038, 702107, 702028, 702732, 702030, 702080, 702081, 702082, 702147, 702402 (PFAS), 702406 (PFAS), 702405 (PFAS)
 N 10 Screw (standard samples)	25	Vials:	702012, 702013
	25	Inserts:	As indicated under Screw N 9 – Inserts
	25	Closures:	702047, 702044, 702046
 N 11 Crimp (standard samples)	26, 27	Vials/	As indicated under GC – Crimp N 11
	26	Inserts/	
	28	Closures:	
 N 11 Snap (standard samples)	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334
	30	Inserts:	As indicated under Screw N 9 – Inserts
	31	Closures:	702064, 702717.2, 702143, 702150, 702151, 702717.2080, 702710, 702710.1, 702095, 702142, 702108, 702710.2080, 702403 (PFAS), 702731, 702063, 702295, 702063.2080
Headspace			HS6, HS40/HS40 Trap/HS40 XL, HS100/101, TurboMatrix™ HS16, HS110/HS110 Trap, Turbomatrix MultiPrep / MultiPrep+, HS2400 M
 N 18 Screw (for Combi PAL® + Turbomatrix MultiPrep (+), TurboMatrix™ HS 16, 40 + 110 and HS2400 M)	37	Vials:	702866, 702866.2, 702826, 702826.0, 702826.2
	37	Closures:	702055, 702827, 702137
 N 20 Crimp (for Combi PAL® / Turbomatrix MultiPrep (+))	38	Vials:	702924, 702263, 702263.0, 702263.2
	41	Closures:	702929, 702834, 702341, 702834, 702834, 702838, 702928.9, 702928
 N 20 Crimp (for TurboMatrix™ HS 16, 40 + 110 and HS2400 M) *** not suitable for TurboMatrix™ 110	38	Vials:	702917***, 70254, 702540
	40, 41,	Closures:	702773, 702829, 70234.8, 702775, 70234.9, 702835, 702094, 702145, 702931.3 + 702804 = 70237.2, 70234.10, 70234
	42		
Purge + Trap (TOC)			AQUATEk 100, Turbomatrix MultiPrep / MultiPrep+
 N 24 Screw	43	Vials:	702023, 702024
	44	Closures:	702058, 702171, 702060, 702062

Autosampler compatibility

Shimadzu®

Main chapter	Pages	Most popular MN products for use on Shimadzu® instruments	
GC			AOC-14, AOC-17, AOC20/AOC-20i/AOC-20s, AOC-30 series, AOC-5000 /AOC5000Plus, AOC-6000
 N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702878
 N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702006, 702007, 702007.0, 702008, 702135, 702335, 702009, 702172, 702010, 702500
	21	Inserts:	702716, 702813, 702077, 702825, 702818, 702819, 702818.1
	22, 23	Closures:	702081, 702732, 702080, 702082, 702147, 702028 (→ all b.m. not for AOC-5000/5000Plus and AOC-6000) 702287.1, 702036, 702037, 702038, 702107, 702402 (PFAS), 702406 (PFAS), 702026 702155, 702156 (for AOC-5000 / AOC5000Plus, AOC-6000)
 N 10 Screw (standard samples)	25	Vials:	702011, 702012, 702013
	25	Inserts:	As indicated under Screw N 9 – Inserts
	25	Closures:	702046, 702048
 N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702892, 702075, 702076, 702888, 702891, 702891.0, 702014, 702134, 702134.1, 702334, 702141
	26	Inserts:	As indicated under Screw N 9 – Inserts
	28	Closures:	702730, 702730.1, 702730.2, 702730.3, (→ all not for AOC-5000/5000Plus and AOC-6000) 70288, 70288.1, 70288.2, 70288.3, 702995, 702995.1, 702995.2, 702995.3, 702879 (for AOC-5000/5000Plus and AOC-6000)
 N 13 Screw (large sample volumes)	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702963 + 702292
HPLC			SIL-2AS, SIL-6A/SIL-6B, SIL-7A, SIL-8A, SIL-9A, LC-2010a/c, SIL-HTa/c, SIL-10A/AF/AP/SIL-10Ai/SILAxL (Prominence), SIL-10ADvp SIL-20A/SIL-20AC (Prominence), SIL-20AXR/SIL-20ACXR (Nexera XR), SIL-30AC/ACMP (Nexera), SIL-40/SIL-40C/SIL-40 XR/SIL-40C XR/SIL-40C XS/SIL-40C X3 (Nexera), SIL-5000
IC*			SIL-20A/AC (HIC-ESP)
 N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702878
 N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702006, 702007, 702007.0, 702008, 702135, 702335, 702009*, 702172*, 702010*, 702500*
	21	Inserts:	702716, 702813, 702077, 702825, 702818, 702819*, 702818.1
	22, 23	Closures:	702287.1, 702036, 702037, 702038, 702107, 702402 (PFAS), 702406 (PFAS), 702405 (PFAS), 702026, 702039, 702040, 702083, 702288.1, 702109, 702027
 N 10 Screw (standard samples)	25	Vials:	702011, 702012, 702013
	25	Inserts:	As indicated under GC – Screw N 9 – Inserts
	25	Closures:	702046, 702047
 N 11 Crimp (standard samples)	26, 27	Vials:	As indicated under GC – Crimp N 11 – Vials
	26	Inserts:	As indicated under GC – Screw N 9 – Inserts
	28	Closures:	70288, 70288.1, 70288.2, 70288.3, 702823
 N 11 Snap (standard samples)	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
	30	Inserts:	As indicated under GC – Screw N 9 – Inserts
	31	Closures:	702710.1, 702710, 702710.2080, 702095, 702142, 702108, 702403 (PFAS), 702717.2, 702064, 702717.2080, 702143, 702150, 702151
 N 13 Screw (large sample volumes)	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702963 + 702292, 702925
 N 8, Shell vials (standard samples)	18	Vials + Plugs:	70202.1 + 702807 (not for Sil-10/Ai/AF/AP, SIL-2AS, SIL-6A/Sil-6B, SIL-7A, SIL-8A, SIL-9A)

Autosampler compatibility

Shimadzu®

Main chapter	Pages	Most popular MN products for use on Shimadzu® instruments	
Headspace		AOC-5000/5000Plus, AOC-6000, HS-10, HS-20/HS-20Trap, HSS-2B, HSS-3A, HT200H	
■ N 18 Screw AOC-5000/5000Plus / AOC-6000, HT200H	37	Vials:	702866, 702866.2, 702826, 702826.0, 702826.2
	37	Closures:	702055, 702827
■ N 20 Crimp (AOC-5000/5000Plus / AOC-6000)	38	Vials:	702924, 702263, 702263.0, 702263.2
	41, 42	Closures:	702929, 702928, 702834, 702341, closure for washer bottle 702924: 70267 + 702144
■ N 20 Crimp (HSS-2B, HSS-3A)	38	Vials:	702918, 702261, 70210.36
	40	Closures:	702093, 702094, 702927, 702835
■ N 20 Crimp (HS-20/HS-20Trap, HT200H, HS-10)	38	Vials:	702924, 702263, 702263.0, 702263.2
	40	Closures:	702093, 702094, 702927, 702835
Purge + Trap (TOC)		ASI-L/ASI-V	
■ N 20 Screw	35	Vials:	702099
	35	Closures:	702181
■ N 24 Screw	43	Vials:	702023, 702024, 702021, 702022
	44	Closures:	702058, 702171, 702060, 702062

Thermo Scientific®

Main chapter	Pages	Most popular MN products for use on Thermo Scientific® instruments	
GC		AI/AS1310, TriPlus™ RSH, TriPlus™ 100 LS, AS1000/AS2000 (Trace GC), AS/AI3000(II) (Trace + Focus GC), AS200, AS800, AI/AS 1610 (Trace GC)	
■ N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702025, 702878
■ N 8 Screw (standard samples)	19	Vials:	70213, 70213.2, 702004, 702893, 702860
	19	Inserts:	702968, 702968.1, 702824, 702005
	20	Closures:	702067, 70245, 702069
■ N 9 Screw (standard samples)	21	Vials:	702282, 702283, 702293, 702284, 702078, 702079, 702006, 702007, 702007.0, 702008, 702135, 702335, 702009, 702172, 702010, 702500
	21	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819, 702825
	22, 23	Closures:	702732, 702030, 702080, 702081, 702082, 702147, 702028, 702287.1, 702287, 702036, 702037, 702038, 702107, 702026, 702402 (for PFAS), 702406 (for PFAS), 702155 (for TriPlus™ RSH)
■ N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702892, 702075, 702076, 702888, 702891, 702891.0, 702014, 702134, 702134.1, 702334
	26	Inserts:	As indicated under Screw N 9 – Inserts
	28	Closures:	70256, 702730, 702730.1, 702730.2, 702730.3, 70288.3, 70288, 70288.2, 70288.1, 702879 (for TriPlus™ RSH)

Autosampler compatibility

Thermo Scientific®

Main chapter	Pages	Most popular MN products for use on Thermo Scientific® instruments	
HPLC		Famos (LC Packings), AS1000/AS3000/AS3500 SpectraSYSTEM, AS100/AS300 SpektraSERIES, Surveyor/Surveyor Plus, Accela/Accela Open AS as well as autosamplers of Dionex® (pls. see page 10)	
■ N 8 Crimp (microsampling)	18	Vials:	70282, 70286
	18	Closures:	70289, 702025, 702878
	19	Vials:	70213, 70213.2, 702004, 702893, 702860
■ N 8 Screw (standard samples)	19	Inserts:	702968, 702968.1, 702824, 702005
	20	Closures:	70245, 702437
	21	Vials:	702282, 702283, 702293, 702284, 702078, 702079, 702006, 702007, 702007.0, 702008, 702135, 702335, 702009, 702172, 702010, 702500
■ N 9 Screw (standard samples)	21	Inserts:	702716, 702813, 702077, 702818, 702818.1, 702819, 702825
	22, 23	Closures:	702287.1, 702287, 702036, 702037, 702038, 702107, 702026, 702402 (PFAS), 702406 (PFAS), 702405 (PFAS), 702288.1, 702288, 702039, 702040, 702083, 702109, 702027
	26, 27	Vials:	As indicated under GC – Crimp N 11 - Vials
■ N 11 Crimp (standard samples)	26	Inserts:	As indicated under Screw N 9 – Inserts
	28	Closures:	70288.3, 70288, 70288.2, 70288.1, 702823
	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
■ N 11 Snap (standard samples)	30	Inserts:	As indicated under GC – Screw N 9 – Inserts
	31	Closures:	702710, 702710.1, 702710.2080, 702095, 702142, 702108, 702403 (PFAS), 702064, 702717.2, 702717.2080, 702150, 702151, 702143
■ N 20 Crimp (large sample volumes)	38	Vials:	702924
	40	Closures:	702093, 702340
Headspace		HS250/500/800/850 (Carlo Erba), HS2000 (Finnigan), TriPlus™ 300 HS, TriPlus™ 500 HS, TriPlus™ RSH, TriPlus™ 100 LS	
■ N 18 Screw (for TriPlus™ RSH/ HS/ SPME/300/500)	37	Vials:	702866, 702866.2, 702826, 702826.0, 702826.2
	37	Closures:	702055, 702827
■ N 20 Crimp (*: for TriPlus™ RSH, TriPlus™ 100LS, HS250/500/800)	38	Vials:	702924*, 702918, 702263*, 702263.0*, 702263.2*, 702261
	40, 41, 42	Closures:	702775, 70234.9, 702773, 702931.3 + 702804 = 70237.2, 702093, 702094, 702340, 70234.10, 70234, 702929*, 702834*, 702341*
Purge + Trap (TOC)		HiPerTOC	
■ N 24 Screw	43	Vials:	702021, 702022, 702023, 702024
	44	Closures:	702058, 702171, 702168, 702060, 702062

Autosampler compatibility

VWR® (Merck® / Hitachi®)

Main chapter	Pages	Most popular MN products for use on VWR® instruments	
HPLC		AS 2000, AS 4000, AS 6000, LaChrom Elite L-2200/L-2200U, LaChrom Elite L-7200/L-7250, Chromaster HPLC 5210/5260/5280, ChromasterUltra RS UHPLC 6270, Primaide™ PM1210 / PM1210TS, 655-A40, L-9100	
■ N 8 Crimp (microsampling) For AS 4000/AS 6000, LaChrom L-7200/L-7250	18	Vials:	70282, 70286
	18	Closures:	70289, 702878
■ N 8 Screw (standard samples)	19	Vials:	70213, 70213.2, 702004, 702893, 702860
	19	Inserts:	702968, 702968.1, 702824, 702005
	20	Closures:	70245, 702437, 702067
■ N 9 Screw (standard samples)	21	Vials:	702282, 702293, 702283, 702078, 702284, 702079, 702007, 702007.0, 702008, 702006, 702135, 702335, 702009, 702172, 702010, 702500
	21	Inserts:	702813, 702077, 702818, 702818.1, 702819, 702716, 702825
	22, 23	Closures:	702287.1, 702287, 702036, 702037, 702038, 702107, 702402 (PFAS), 702406 (PFAS), 702405 (PFAS), 702288.1, 702288, 702039, 702040, 702083, 702109, 702026, 702027
■ N 11 Crimp (standard samples)	26, 27	Vials:	70201HP, 70201HP.2, 702885, 702075, 702892, 702076, 702888, 702891, 702891.0, 702014, 702134, 702134.1, 702334
	26	Inserts:	As indicated under Screw N 9 – Inserts
	28	Closures:	70288, 70288.1, 70288.2, 70288.3, 702823
■ N 11 Snap (standard samples)	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
	30	Inserts:	As indicated under Screw N 9 – Inserts
	31	Closures:	702710.1, 702710, 702710.2080, 702095, 702142, 702108, 702403 (PFAS), 702717.2, 702064, 702717.2080, 702143, 702150, 702151, 702718.1, 702718, 702718.2080, 702063, 702731, 702063.2080, 702295
■ N 8 Shell vials (standard samples) Only for Chromaster and LaChrom Elite L-2200/ L-2200U	18	Vials + plugs:	70202.1 + 702807
■ N 13 Screw (large sample volumes) Not for AS 2000/AS 4000/ AS 6000	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702925, 702963 + 70260

Waters®

Main chapter	Pages	Most popular MN products for use on Waters® instruments	
HPLC		WISP, 717 (Plus), GPC 2000, 2707/2707 LSV, Alliance® 2690/2695/e2695, Alliance® 2790/2795, ACQUITY® UPLC/UHPLC Systems, Arc HPLC System, 2777 (CTC PAL®)	
■ N 9 Screw (standard samples) Not for WISP, 717 (Plus), GPC 2000	21	Vials:	702282, 702293, 702283, 702284, 702078, 702079, 702007, 702007.0, 702008, 702135, 702335, 702006 (not for 2690/2695), 702009, 702172, 702010, 702500
	21	Inserts:	702818, 702818.1, 702819
	22, 23	Closures:	702026, 702027, 702287.1, 702287, 702036, 702037, 702038, 702402 (PFAS), 702406 (PFAS), 702405 (PFAS), 702288.1, 702288, 702039, 702040, 702083
■ N 10 Screw (standard samples) For Alliance 2690/2695, Alliance 2790/2795	25	Vials:	702011, 702012, 702013
	25	Inserts:	702818, 702818.1, 702819
	25	Closures:	702046, 702047
■ N 11 Crimp (standard samples) For Alliance 2690/2695	26	Vials:	70201HP, 70201HP.2, 702885, 702892, 702075, 702076
	26	Inserts:	702818, 702818.1, 702819
	28	Closures:	70288, 70288.1, 70288.2, 70288.3, 702995, 702995.1, 702995.2, 702995.3
■ N 11 Snap (standard samples) For Alliance 2690/2695, Alliance 2790/2795	30	Vials:	702714, 702713, 702712, 702709, 702708, 702134, 702134.1, 702334, 702809, 702173, 702174
	30	Inserts:	702818, 702818.1, 702819
	31	Closures:	702710.1, 702095, 702142, 702108, 702403 (PFAS), 702717.2, 702143, 702150, 702151
■ N 8 Shell Vials (standard samples) Only for WISP, 717 (Plus)	18	Vials + Plugs:	70202.1 + 702807
■ N 13 Screw (large sample volumes) Only for WISP, 717 (Plus), GPC 2000	33	Vials:	702962, 702973
	33	Insert:	702972 + spring 702974
	33	Closures:	702926, 702925, 702963 + 70260



Finding instead of searching:

Translation of cross-references with the VialFinder



Are you in need to change your current vials and caps supplier due to quality, price or delivery issues? Make an easy changeover to vials and caps from MACHEREY-NAGEL by using our VialFinder – a database driven tool on our website that can help you with the translation of cross-references of vials and caps products from suppliers from all over the world. This powerful and easy-to-use tool has been updated, so that you now find cross-references from almost 60 suppliers worldwide (instrument companies and consumables suppliers).



www.mn-net.com/vialfinder

Shell vials N 8

Fire-polished neck
N 8To be used with
N 8 PE plugs

70202.1

702807

1 mL

8.2 × 40 mm

clear

flat bottom

100 / PP box

100 / PE Bag

Crimp neck vials N 8

For micro sampling



70286

70282

0.2 mL

0.3 mL

5.5 × 31.5 mm

5.5 × 31.5 mm

clear

clear

conical

round bottom

100 / PE bag

100 / PE bag

Crimping tools N 8

REF	Type of crimping tool	Pack of
735126	Manual crimper (standard) for 8 mm aluminum crimp caps	1 / box
735408	Manual decapper (standard) for 8 mm aluminum crimp caps	1 / box
735208	Manual ergonomic crimper for 8 mm aluminum crimp caps	1 / box



Ready assembled crimp closures N 8 and plain crimp caps N 8

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702025		N 8 aluminum crimp cap, silver, center hole	Red Rubber / FEP colorless	45° shore A	1.0 mm	100 / PE bag
70289		N 8 aluminum crimp cap, silver, center hole	Silicone white / PTFE red	45° shore A	1.0 mm	100 / PE bag
702878		N 8 aluminum crimp cap, silver, center hole	PTFE red / Silicone white / PTFE red	50° shore A	1.0 mm	100 / PE bag
702800		N 8 aluminum crimp cap, silver, center hole (no liner)	–	–	–	100 / PE bag



Inserts and micro-vials for sample volumes < 1.0 mL

Samples with a lower volume than 1.0 mL have to be delivered to the autosampler in appropriate micro sample containers, in order to offer the needle ideal conditions for sample removal. During selection of the most suitable product, safety, handling and price are important factors. In order to give you an overview on the different options on micro-inserts and micro-vials, please watch our video "Micro-sampling in chromatography" under: www.mn-net.com/Chromatography/Vials-and-caps/Inserts/



Screw neck vials N 8

Small opening, 8-425 thread

				
70213	70213.2	702004	702893	702860
1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.1 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm
clear	amber	clear, label + scale	amber, label + scale	clear
flat bottom	flat bottom	flat bottom	flat bottom	conical
100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PE bag

Inserts for screw neck vials N 8 (small opening)

Inserts for small opening vials

			
702968*	702968.1*	702824	702005
0.1 mL	0.15 mL	0.1 mL	0.25 mL
5 x 31 mm	5 x 31 mm	5 x 29 mm	5 x 31 mm
clear	clear	clear	clear
conical, 15 mm tip	conical, 9 mm tip	with plastic spring	flat bottom
100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag

Container for screw neck vials N 8

REF	Description	Pack of
702514	81 position container blue for screw neck vials N 8, 130 x 130 x 45 mm, coded, with transparent lid (suitable for freezers)	1

Vial rack for screw neck vials N 8

REF	Description	Pack of
702502	50 position polypropylene vial rack blue for all vials 11.6 x 32 mm with flat bottom Dimensions: 190 x 100 x 22 mm, stackable	1



N 8 Screw neck vials, inserts and caps N 8

Ready assembled screw closures N 8 and plain screw caps N 8

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702067		N 8 PP screw cap, black, center hole	Red Rubber / FEP colorless	45° shore A	1.3 mm	100 / PE bag
702068		as above, but with closed top				
70245		N 8 PP screw cap, black, center hole	Silicone white / PTFE red	45° shore A	1.3 mm	100 / PE bag
702066		as above, but with closed top				
702437		N 8 PP screw cap, black, center hole	Silicone white / PTFE blue, slit	55° shore A	1.0 mm	100 / PE bag
702069		N 8 PP screw cap, black, center hole	PTFE red / Silicone white / PTFE red	50° shore A	1.0 mm	100 / PE bag
70249		N 8 PP screw cap, black, center hole (no liner)	–	–	–	100 / PE bag
70250		N 8 PP screw cap, black, closed top (no liner)	–	–	–	100 / PE bag

N 8 septa for screw caps N 8

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
70261		- (septa only)	PTFE virginal, white	53° shore D	0.25 mm	100 / PE bag
702070		- (septa only)	Red Rubber / FEP colorless	45° shore A	1.3 mm	100 / PE bag
70248		- (septa only)	Silicone white / PTFE red	45° shore A	1.3 mm	100 / PE bag
702481		- (septa only)	Silicone white / PTFE blue, slit	55° shore A	1.0 mm	100 / PE bag



Storage of samples in the fridge or in the freezer

Useful tips for sample handling

Generally sample vials should be stored in a vial container when being placed in the fridge or in the freezer, in order to avoid any condensations on the cap/septa surface that may go along with contaminations in the penetration area of the septa in the center hole. When filling the vial you have to consider the expansion rate of your sample to prevent breakage of the vial. Furthermore it is important to defreeze the sample at a later point in time very slowly (no sudden defreezing with hot water for example). With screw closures you may have to check, if restoring forces have been activated during the defreezing process and if you may have to tighten the screw closure. The choice of the correct closure (septum) depends on the storage temperature.



Screw neck vials and inserts N 9

Wide opening, short thread									
Polypropylene					Polypropylene				
									
702282	702293	702283	702284	702500	702813	702716	702818	702819	702825
		silanized 702078	silanized 702079		silanized 702077		silanized 702818.1		
1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL	0.2 mL	0.25 mL	0.1 mL	0.1 mL	0.3 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	6 x 31 mm	6 x 31 mm	5.7 x 29 mm	5.7 x 29 mm	6 x 31 mm
clear	amber	clear, label + scale	amber, label + scale	transparent, filling lines	clear	clear	clear	transparent	clear
flat bottom	flat bottom	flat bottom	flat bottom	flat bottom	conical, 15 mm tip	conical, 12 mm tip	with plastic spring	with integrated spring	flat bottom
100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag

Micro-vials N 9 for small sample volumes									
				PP / Glass insert	PP / Glass insert		Polypropylene	Polypropylene	Polypropylene
									
702006	702007	702007.0	702008	702135*	702335*	702088	702009	702172	702010
1.1 mL	0.2 mL	0.2 mL	0.2 mL	0.2 mL	0.2 mL	1.1 mL	0.3 mL	0.3 mL	0.7 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm
clear	clear	clear	amber	clear	amber	clear	transparent	amber	transparent
flat bottom, 15 µL funnel in solid glass bottom	flat bottom, with integrated 0.2 mL insert	flat bottom, with integra- ted 0.2 mL insert, label	flat bottom, with integra- ted 0.2 mL insert	flat bottom, with integra- ted 0.2 mL glass insert	flat bottom, with integra- ted 0.2 mL glass insert	conical with round pedestal glass plate	with inner cone	with inner cone	with round bottom insert
100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PE bag	100 / PE bag	100 / PE bag

* upon request also available with an integrated silanized insert

Pre-assembled vial-insert combinations with screw neck N 9

Vial description	Insert description	Pack of	REF
Vial 702282: 1.5 mL, clear, flat bottom	with pre-assembled micro-insert 702813: 0.2 mL, conical, 15 mm tip	100 / PP box	702177
Vial 702283: 1.5 mL, clear, flat bottom, label + scale	with pre-assembled micro-insert 702813: 0.2 mL, conical, 15 mm tip	100 / PP box	702178
Vial 702284: 1.5 mL, amber, flat bottom, label + scale	with pre-assembled micro-insert 702813: 0.2 mL, conical, 15 mm tip	100 / PP box	702179
Further pre-assembled vial-insert combinations on request			

Container for screw neck vials N 9

REF	Description	Pack of
702514	81 position container blue for screw neck vials N 9, 130 x 130 x 45 mm, coded, with transparent lid (suitable for freezers)	1

Vial rack for screw neck vials N 9

REF	Description	Pack of
702502	50 position polypropylene vial rack blue for all vials 11.6 x 32 mm with flat bottom Dimensions: 190 x 100 x 22 mm, stackable	1



Ready assembled screw closures, plain screw caps and single septa N 9

Cap description (Scale 1:1.6)	Septum description	Hardness	Thickness	Pack of
N 9 PP screw cap, color as indicated, center hole   702029 702031	closed top PTFE virginal, white 	53° shore D	0.25 mm	100 / PE bag
 702032				
N 9 PP screw cap, blue, center hole   702402 702405	Silicone white / Polyimide orange, (fluorine-free, for PFAS analysis) as above, however, with slit  	50° shore A	1.0 mm	100 / PE bag
N 9 PP screw cap, blue, center hole  702406	Silicone white / Polypropylene colorless (fluorine-free, for PFAS analysis) 	45° shore A	1.0 mm	100 / PE bag
N 9 PP screw cap, color as indicated, center hole        702030 702732 702080 702081 702082 702147 702033	closed top Red Rubber / FEP colorless 	45° shore A	1.0 mm	100 / PE bag
N 9 PP screw cap, color as indicated, center hole        702287 702287.1 702036 702037 702038 702107 702034	closed top Silicone white / PTFE red 	45° shore A	1.0 mm	100 / PE bag
N 9 PP screw cap, color as indicated, center hole       702288 702288.1 702039 702040 702083 702109	Silicone white / PTFE blue, slit 	55° shore A	1.0 mm	100 / PE bag
N 9 PP screw cap, color as indicated, center hole       702286 702035 702158 702084 702085 702159	PTFE red / Silicone white / PTFE red 	50° shore A	1.0 mm	100 / PE bag
N 9 PP screw cap, color as indicated, center hole (no liner)       702160 702161 702162 702163 702164 702165	–	–	–	100 / PE bag

Ready assembled magnetic screw closures N 9

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702155		N 9 magnetic screw cap, silver, center hole	Silicone white / PTFE red	45° shore A	1.0 mm	100 / PE bag
702156		N 9 magnetic screw cap, silver, center hole	Silicone white / PTFE blue, slit	55° shore A	1.0 mm	100 / PE bag

N 9 septa for screw caps N 9

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702043		- (septa only)	PTFE virginal, white	53° shore D	0.25 mm	100 / PE bag
702041		- (septa only)	Red Rubber / FEP colorless	45° shore A	1.0 mm	100 / PE bag
702042		- (septa only)	Silicone white / PTFE red	45° shore A	1.0 mm	100 / PE bag
702148		- (septa only)	Silicone white / PTFE blue, slit	55° shore A	1.0 mm	100 / PE bag

Bonded screw closures N 9

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702028		N 9 PP bonded screw cap, blue, center hole	Red Rubber / TEF colorless	65° shore A	1.0 mm	100 / PE bag
702026		N 9 PP bonded screw cap, blue, center hole	Silicone beige / PTFE white	45° shore A	1.3 mm	100 / PE bag
702027		N 9 PP bonded screw cap, blue, center hole	Silicone beige / PTFE white, slit	45° shore A	1.3 mm	100 / PE bag

(Septum firmly connected with the cap; cannot be removed)



Vials and caps for PFAS analysis

Per- and polyfluoroalkyl substances (PFAS) are used for many daily applications, e.g. food packaging and textile coating, because of their non-sticky properties. However, most PFAS are dangerous substances that need to be monitored globally. Due to their increasing appearance in the environment, new product solutions for their analysis are of great demand.

When you are doing PFAS analysis, it is crucial to select the right vials and closures for this application. Adsorption effects of glass as well as possible contaminations of the sample by particles from the septa, especially from the PTFE lamination, may put your analysis results at risk.

Polypropylene vials are best suited for PFAS analysis, since the adsorption effects are the lowest and the signal strength of the analytes is therefore the highest. Fluorine-free septa, i. e. septa with a polyimide or a polypropylene coating rather than a PTFE coating, are recommended to eliminate any migration of fluorine into the sample.



N9 Screw neck vials, inserts and caps N 9



Pre-sealed vial-closure combinations with screw neck N 9

Vial description	Closure description	Pack of	REF
Pre-sealed vials 702282 : 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening	pre-screwed with 702732 : N 9 PP screw cap, blue, center hole, Red Rubber / FEP colorless, 45° shore A, 1.0 mm	100	702857
Pre-sealed vials 702283 : 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening, label + scale	pre-screwed with 702732 : N 9 PP screw cap, blue, center hole, Red Rubber / FEP colorless, 45° shore A, 1.0 mm	100	702858
Pre-sealed vials 702282 : 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening	pre-screwed with 702287.1 : N 9 PP screw cap, blue, center hole, Silicone white / PTFE red, 45° shore A, 1.0 mm	100	702874
Pre-sealed vials 702283 : 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening, label + scale	pre-screwed with 702288.1 : N 9 PP screw cap, blue, center hole, Silicone white / PTFE blue, slit, 55° shore A, 1.0 mm	100	702863
Pre-sealed vials 702284 : 1.5 mL screw neck vial N 9, 11.6 x 32 mm, amber, flat bottom, wide opening, label + scale	pre-screwed with 702288.1 : N 9 PP screw cap, blue, center hole, Silicone white / PTFE blue, slit, 55° shore A, 1.0 mm	100	702873
Pre-sealed vials 702283 : 1.5 mL screw neck vial N 9, 11.6 x 32 mm, clear, flat bottom, wide opening, label + scale	pre-screwed with 702026 : N 9 PP bonded screw cap, blue, center hole, Silicone beige / PTFE white, 45° shore A, 1.3 mm	100	702864
Other pre-sealed vial-closure combinations on request			



Special vials for special applications

Silanized glass vials / Plastic vials / Plastic vials with glass insert

▪ Silanized glass vials

Silanized glass vials have a de-activated inner glass surface, in order to reduce adsorption of polar substances. Therefore they are often used for the analysis of proteins, phenols and amino acids, which would – without any silanization of the glass surface – react with the OH-groups of the glass and thus would stick to the normally polar glass surface. It is also recommendable to use silanized vials respectively inserts for pH-sensitive and aqueous samples.

▪ Plastic vials

For some applications glass vials are not suitable due to their composition and their chemical properties. Amongst these are heavy metal analysis, water and protein analysis, atomic absorption, capillary electrophoresis (CE) and ion chromatography (IC). For all these cases high purity polypropylene vials with 1.5 mL, 0.3 mL and 0.7 mL volume in transparent and amber are available.

▪ Plastic vials with integrated glass insert

In comparison to the glass-in-glass products, glass-in-plastic systems are very robust, as the glass insert is well protected by the polypropylene outer shell. The tip of the micro-insert is centered by 100 per cent in an outlet at the bottom. The insert sits firmly in the protective PP round bottom shell and thus can easily be filled. Another advantage of these systems is their excellent tightness, as the glass insert always constantly exceeds the rim of the plastic outer vial by 0.1 mm granting a firm sealing of the sample in the insert. Upon request also a silanized insert can be integrated into the plastic shell. The high transparent polypropylene enables a good view on the filling level.



Screw neck vials and inserts N 10

Wide opening, 10-425 thread

Polypropylene							
							
702011	702012	702013	702813	702716	702818	702819	702825
1.5 mL	1.5 mL	1.5 mL	silanized 702077		silanized 702818.1		
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	0.2 mL	0.25 mL	0.1 mL	0.1 mL	0.3 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	6 x 31 mm	6 x 31 mm	5.7 x 29 mm	5.7 x 29 mm	6 x 31 mm
clear	clear, label + scale	amber, label + scale	clear	clear	clear	transparent	clear
flat bottom	flat bottom	flat bottom	conical, 15 mm tip	conical, 12 mm tip	with plastic spring	with integrated spring	flat bottom
100 / PP box	100 / PP box	100 / PP box	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag

Container for screw neck vials N 10

REF	Description	Pack of
702514	81 position container blue for screw neck vials N 10, 130 x 130 x 45 mm, coded, with transparent lid (suitable for freezers)	1

Vial rack for screw neck vials N 10

REF	Description	Pack of
702502	50 position polypropylene vial rack blue for all vials 11.6 x 32 mm with flat bottom Dimensions: 190 x 100 x 22 mm, stackable	1



Screw closures N 10 and plain screw caps N 10

REF	Scale (1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702044		N 10 PP bonded* screw cap, black, center hole	Red Rubber / TEF colorless	65° shore A	1.0 mm	100 / PE bag
702046		N 10 PP bonded* screw cap, black, center hole	Silicone white / PTFE red	45° shore A	1.0 mm	100 / PE bag
702047		N 10 PP bonded* screw cap, black, center hole	Silicone white / PTFE blue, slit	45° shore A	1.5 mm	100 / PE bag
702048		N 10 PP screw cap, black, center hole (ready assembled)	PTFE red / Silicone white / PTFE red	50° shore A	1.0 mm	100 / PE bag
702049		N 10 PP screw cap, black, center hole (no liner)	-	-	-	100 / PE bag

* Septum firmly connected with the cap; cannot be removed

N 11 Crimp neck vials, inserts and caps N 11

Crimp neck vials N 11

Small opening vials		Wide opening vials			
					
70201CG	70214CG	70201HP	70201HP.2	702885 silanized 702075	702892 silanized 702076
1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm
clear	amber	clear	amber	clear, label + scale	amber, label + scale
flat bottom	flat bottom	flat bottom	flat bottom	flat bottom	flat bottom
100/PP box	100/PP box	100/PP box	100/PP box	100/PP box	100/PP box

Inserts for crimp neck vials N 11

Inserts for small opening vials				Inserts for wide opening vials				
								
702968	702968.1	702824	702005	702813 silanized 702077	702716	702818 silanized 702818.1	702819	702825
0.1 mL	0.15 mL	0.1 mL	0.25 mL	0.2 mL	0.25 mL	0.1 mL	0.1 mL	0.3 mL
5 x 31 mm	5 x 31 mm	5 x 29 mm	5 x 31 mm	6 x 31 mm	6 x 31 mm	5.7 x 29 mm	5.7 x 29 mm	6 x 31 mm
clear	clear	clear	clear	clear	clear	clear	transparent	clear
conical, 15 mm tip	conical, 9 mm tip	with plastic spring	flat bottom	conical, 15 mm tip	conical, 12 mm tip	with plastic spring	with integrated spring	flat bottom
100/PE bag	100/PE bag	100/PE bag	100/PE bag	100/PE bag	100/PE bag	100/PE bag	100/PE bag	100/PE bag

Container for crimp neck/snap ring vials N 11

REF	Description	Pack of
702514	81 position container blue for crimp neck/snap ring vials N 11, 130 x 130 x 45 mm, coded, with transparent lid (suitable for freezers)	1

Vial rack for crimp neck/snap ring vials N 11

REF	Description	Pack of
702502	50 position polypropylene vial rack blue for all vials 11.6 x 32 mm with flat bottom Dimensions: 190 x 100 x 22 mm, stackable	1



Crimp neck vials N 11

Micro-vials N 11 for small sample volumes

							PP / Glass insert	PP / Glass insert
								
702888	702141	702015	702016	702891	702891.0	702014	702134*	702334*
1.1 mL	1.1 mL	1.1 mL	1.1 mL	0.2 mL	0.2 mL	0.2 mL	0.2 mL	0.2 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm
clear	clear	clear	amber	clear	clear	amber	clear	amber
flat bottom, 15 µL funnel in solid glass bottom	conical	conical with round pedestal glass plate	conical with round pedestal glass plate	flat bottom, with integrated 0.2 mL insert	flat bottom, with integrated 0.2 mL insert, label	flat bottom, with integrated 0.2 mL insert	flat bottom, with integrated 0.2 mL glass insert	flat bottom, with integrated 0.2 mL glass insert
100 / PP box	100 / PE bag	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box

* upon request also available with integrated silanized insert



Crimping tools N 11

REF	Type of crimping tool	Pack of
735111	Manual crimper (standard), height adjustable, for 11 mm aluminum crimp caps	1 / box
735911	Manual decapper (standard) for 11 mm aluminum crimp caps	1 / box
735211	Ergonomic crimper for 11 mm aluminum crimp caps	1 / box
735311	Ergonomic decapper for 11 mm aluminum crimp caps	1 / box
735511	Electronic crimper for 11 mm aluminum crimp caps (battery-powered)	1 / box
735611	Electronic decapper for 11 mm aluminum crimp caps (battery-powered)	1 / box
735700	Electronic high power crimping tool with power supply	1 / box
735711	Crimping head for 11 mm crimp caps (aluminum, magnetic)	1 / box
735811	Decapping head for 11 mm crimp caps (aluminum, magnetic)	1 / box



N 11 Crimp neck vials, inserts and caps N 11

Ready assembled crimp closures and plain crimp caps N 11

Cap description (Scale 1:1.6)	Septum description	Hardness	Thickness	Pack of
N 11 aluminum crimp cap, color as indicated, center hole  70246 70246.1 70246.2 70246.3	Natural rubber orange / ETFE colorless 	45° shore A	1.0 mm	100 / PE bag
N 11 aluminum crimp cap, color as indicated, center hole  70256 70231.1 70231.2	Natural rubber / Butyl red-orange / TEF colorless 	45° shore A	1.0 mm	100 / PE bag
N 11 aluminum crimp cap, color as indicated, center hole  702730 702730.1 702730.2 702730.3	Red Rubber / FEP colorless 	45° shore A	1.0 mm	100 / PE bag
N 11 aluminum crimp cap, color as indicated, center hole  70288 70288.1 70288.2 70288.3 702879 (magnetic, 1.0 mm)	Silicone white / PTFE red 	45° shore A	1.3 mm	100 / PE bag
N 11 aluminum crimp cap, color as indicated, center hole  702995 702995.1 702995.2 702995.3	PTFE red / Silicone white / PTFE red 	50° shore A	1.0 mm	100 / PE bag
N 11 aluminum crimp cap, silver, center hole  702823*	Silicone white / PTFE blue, cross-slit 	55° shore A	1.5 mm	100 / PE bag
N 11 aluminum crimp cap, silver, center hole  702736*	Red Rubber / FEP colorless 	45° shore A	1.3 mm	100 / PE bag
N 11 aluminum crimp cap, silver, center hole  702001	Natural rubber red-orange / TEF colorless 	60° shore A	1.0 mm	100 / PE bag
N 11 aluminum crimp cap, silver, center hole  70231	Natural rubber / Butyl red-orange / TEF colorless 	45° shore A	1.3 mm	100 / PE bag
N 11 aluminum crimp cap, silver, center hole  70239	PTFE gray / Butyl beige / PTFE gray 	55° shore A	1.3 mm	100 / PE bag
N 11 aluminum crimp cap, silver, center hole  70284	PTFE virginal, white 	53° shore D	0.25 mm	100 / PE bag
N 11 aluminum crimp cap, silver, roll groove, center hole  702175	O-ring and aluminum septa TPF (Total Phthalate Free) 		0.1 mm	100 / glass jar

* upon request also available with a green, red or a blue crimp cap

Ready assembled crimp closures and plain crimp caps N 11

Cap description (Scale 1:1.6)	Septum description	Hardness	Thickness	Pack of
N 11 aluminum crimp cap, silver, center hole (no liner)	–	–	–	100/PE bag
				
702801				
N 11 PE cap, transparent, closed top, thin piercing area	–			100/PE bag
				
702401				

* upon request also available with a green, red or a blue crimp cap

Pre-sealed vial-closure combinations with crimp neck N 11

Vial description	Closure description	Pack of	REF
Pre-sealed vials 70201HP : 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, clear, flat bottom, wide opening	crimped with 70256 : N 11 aluminum crimp cap, silver, center hole, Natural rubber / Butyl red-orange / TEF colorless, 45° shore A, 1.0 mm	100	702101HP
Pre-sealed vials 702892 : 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, amber, flat bottom, wide opening, label and scale	crimped with 70256 : N 11 aluminum crimp cap, silver, center hole, Natural rubber / Butyl red-orange / TEF colorless, 45° shore A, 1.0 mm	100	702859
Pre-sealed vials 70201HP : 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, clear, flat bottom, wide opening	crimped with 702995 : N 11 aluminum crimp cap, silver, center hole, PTFE red / Silicone white / PTFE red, 50° shore A, 1.0 mm	100	702867
Pre-sealed vials 70201CG : 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, clear, flat bottom, small opening	crimped with 70231 : N 11 aluminum crimp cap, silver, center hole, Natural rubber / Butyl red-orange / TEF colorless, 45° shore A, 1.3 mm	100	702882
Pre-sealed vials 70201HP : 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, clear, flat bottom, wide opening	crimped with 702823 : N 11 aluminum crimp cap, silver, center hole, Silicone white / PTFE blue, cross-slit, 55° shore A, 1.5 mm	100	702887
Pre-sealed vials 702892 : 1.5 mL crimp neck vial N 11, 11.6 x 32 mm, amber, flat bottom, wide opening, label + scale	crimped with 702823 : N 11 aluminum crimp cap, silver, center hole, Silicone white / PTFE blue, cross-slit, 55° shore A, 1.5 mm	100	702895
Other pre-sealed vial-closure combinations on request			



Useful tips for optimal crimping

can be found on our website under www.mn-net.com/media/pdf/Poster-Optimal-Crimping-EN.pdf.

As a decision-making aid when choosing the most suitable crimping tool for your personal needs, use our specialized brochure “Crimping tools”, which you can find in the chromatography download area at www.mn-net.com/chromatography-downloads on our website.

To find out about the different operation modes of the various crimping tool types, please click on the tab “Learn more” on the product detail page of each crimping tool on our website.



N 11 Snap ring vials, inserts and caps N 11

Snap ring vials and inserts N 11

Wide opening, vials can alternatively be crimped as well (see page 28)

Polypropylene							
							
702714	702713	702712	702813 silanized 702077	702716	702818 silanized 702818.1	702819	702825
1.5 mL	1.5 mL	1.5 mL	0.2 mL	0.25 mL	0.1 mL	0.1 mL	0.3 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	6 x 31 mm	6 x 31 mm	5.7 x 29 mm	5.7 x 29 mm	6 x 31 mm
clear	clear, label + scale	amber, label + scale	clear	clear	clear	transparent	clear
flat bottom	flat bottom	flat bottom	conical, 15 mm tip	conical, 12 mm tip	with plastic spring	with integrated spring	flat bottom
100 / PP box	100 / PP box	100 / PP box	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag

Micro-vials N 11 for small sample volumes


	PP / Glass insert	PP / Glass insert	Polypropylene	Polypropylene	Polypropylene
					
702709/702708	702134*	702334*	702809	702173	702174
0.2 mL	0.2 mL	0.2 mL	0.3 mL	0.3 mL	0.7 mL
11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm	11.6 x 32 mm
clear/amber	clear	amber	transparent	amber	transparent
flat bottom, with integrated 0.2 mL insert	flat bottom, with integrated 0.2 mL glass insert	flat bottom, with integrated 0.2 mL glass insert	with inner cone	with inner cone	with round bottom insert
100 / PP box	100 / PP box	100 / PP box	100 / PE bag	100 / PE bag	100 / PE bag

* upon request also available with an integrated silanized insert

Pre-assembled vial-insert combinations with snap ring N 11

Vial description	Insert description	Pack of	REF
Vial 702714: 1.5 mL, clear, flat bottom	with pre-assembled micro-insert 702813: 0.2 mL, conical, 15 mm tip	100 / PP box	702170
Vial 702713: 1.5 mL, clear, flat bottom, label + scale	with pre-assembled micro-insert 702813: 0.2 mL, conical, 15 mm tip	100 / PP box	702176
Further pre-assembled vial-insert combinations on request			

Ready assembled snap ring closures N 11 (hard cap)

Cap description (Scale 1:1.6)	Septum description	Hardness	Thickness	Pack of
N 11 PE snap ring cap, color as indicated, center hole  702731 702063 702295	Red Rubber / FEP colorless 	45° shore A	1.0 mm	100 / PE bag
N 11 PE snap ring cap, color as indicated, center hole  702710 702710.1 702095 702142 702108	Silicone white / PTFE red 	45° shore A	1.0 mm	100 / PE bag
N 11 PE snap ring cap, color as indicated, center hole  702064 702717.2 702143 702150 702151	Silicone white / PTFE blue, cross-slit 	55° shore A	1.0 mm	100 / PE bag
N 11 PE snap ring cap, color as indicated, center hole  702718 702718.1	PTFE red / Silicone white / PTFE red 	50° shore A	1.0 mm	100 / PE bag
N 11 PE cap, transparent, closed top, thin piercing area  702401	–	–	–	100 / PE bag

Ready assembled snap ring closures N 11 (soft cap)

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702063.2080	 	N 11 PE snap ring cap, soft, light blue, center hole	Red Rubber / FEP colorless	45° shore A	1.0 mm	100 / PE bag
702403	 	N 11 PE snap ring cap, soft, light blue, center hole	Silicone white / Polyimide orange (fluorine-free, for PFAS analysis)	50° shore A	1.0 mm	100 / PE bag
702710.2080	 	N 11 PE snap ring cap, soft, light blue, center hole	Silicone white / PTFE red	45° shore A	1.0 mm	100 / PE bag
702717.2080	 	N 11 PE snap ring cap, soft, light blue, center hole	Silicone white / PTFE blue, cross-slit	55° shore A	1.0 mm	100 / PE bag
702718.2080	 	N 11 PE snap ring cap, soft, light blue, center hole	PTFE red / Silicone white / PTFE red	50° shore A	1.0 mm	100 / PE bag



Inserts and micro-vials for sample volumes < 1.0 mL

Samples with a lower volume than 1.0 mL have to be delivered to the autosampler in appropriate micro sample containers, in order to offer the needle ideal conditions for sample removal. During selection of the most suitable product, safety, handling and price are important factors. In order to give you an overview on the different options on micro-inserts and micro-vials, please watch our video “Micro-sampling in chromatography” under: www.mn-net.com/Chromatography/Vials-and-caps/Inserts/



N 13 Crimp neck vials and caps N 13

Crimp neck vial N 13

70203

2 mL

13.75 x 35 mm

clear, flat bottom

100 / carton



Crimping tools N 13

REF	Type of crimping tool	Pack of
735113	Manual crimper (standard), height adjustable, for 13 mm aluminum crimp caps	1 / box
735133	Manual crimper (standard), height adjustable, for 13 mm Flip Top / Flip Off caps	1 / box
735913	Manual decapper (standard) for 13 mm aluminum crimp caps	1 / box
735513	Electronic crimper for 13 mm aluminum crimp caps (battery-powered)	1 / box
735533	Electronic crimper for 13 mm Flip Top / Flip Off caps (battery-powered)	1 / box
735613	Electronic decapper for 13 mm aluminum crimp caps (battery-powered)	1 / box
735700	Electronic high power crimping tool with power supply	1 / box
735713	Crimping head for 13 mm crimp caps	1 / box
735733	Crimping head for 13 mm Flip Top / Flip Off caps	1 / box
735813	Decapping head for 13 mm crimp caps	1 / box



Ready assembled crimp closures N 13 and plain crimp caps N 13


REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
70257		N 13 aluminum crimp cap, silver, center hole	Butyl dark gray / PTFE gray (only centrally laminated with PTFE, typically called Pharma-Fix)	50° shore A	2 mm	100 / PE bag
70232		N 13 aluminum center tear off cap, gold	Butyl dark gray / PTFE gray (only centrally laminated with PTFE, typically called Pharma-Fix)	50° shore A	2 mm	100 / PE bag
702802		N 13 aluminum crimp cap, silver, center hole (no liner)	-	-	-	100 / PE bag
702803		N 13 aluminum center tear off cap, gold (no liner)	-	-	-	100 / PE bag

Container for crimp and screw neck vials N 13

REF	Description	Pack of
702515	49 position container blue for crimp and screw neck vials N 13, 130 x 130 x 50 mm, with transparent lid (suitable for freezers)	1

Vial rack for crimp and screw neck vials N 13

REF	Description	Pack of
702504	50 position polypropylene vial rack blue for vials with a diameter of 15 mm max. and flat bottom Dimensions: 240 x 120 x 28 mm, stackable	1



Screw neck vials and insert N 13

13-425 Thread



702962	702973	702972	702974
4 mL	4 mL	0.3 mL	
14.75 x 45 mm	14.75 x 45 mm	6 x 40 mm	
clear	amber	clear	
flat bottom	flat bottom	conical, metal spring required	metal spring for 702972
100 / PP box	100 / PP box	100 / PE bag	100 / PE bag

Ready assembled screw closures N 13, plain screw caps N 13, single septa N 12

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702113		N 13 PP bonded* screw cap, white, closed top	Silicone white / PTFE blue	55° shore A	1.3 mm	100 / PE bag
702050		N 13 PP screw cap, black, center hole	Red Rubber / FEP colorless	40° shore A	1.5 mm	100 / PE bag
702051		as above, but with closed top				
702926		N 13 PP screw cap, black, center hole	Silicone white / PTFE red	40° shore A	1.3 mm	100 / PE bag
702052		as above, but with closed top				
702925		N 13 PP screw cap, black, center hole	Silicone white / PTFE blue, cross-slit	50° shore A	1.5 mm	100 / PE bag
702963		N 13 PP screw cap, black, center hole (no liner)	-	-	-	100 / PE bag
702966		as above, but with closed top (no liner)				
70260		- (septa only)	PTFE virginal, white	53° shore D	0.25 mm	100 / PE bag
702053		- (septa only)	Red Rubber / FEP colorless	40° shore A	1.5 mm	100 / PE bag
702292		- (septa only)	Silicone white / PTFE red	40° shore A	1.3 mm	100 / PE bag

* Septum firmly connected with the cap; cannot be removed

Screw neck vials and closed top screw closures for storage of liquid samples

N 8				N 9				N 13		
Small opening, 8-425 thread				Wide opening, short thread				13-425 thread		
70213	70213.2	702004	702893	702282	702293	702283 silanized 702078	702284 silanized 702079	702962	702973	
1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL	1.5 mL	4 mL	4 mL	
11.6 x 32 mm				11.6 x 32 mm				14.75 x 45 mm	14.75 x 45 mm	
clear, flat bottom	amber, flat bottom	clear, flat bottom, label + scale	amber, flat bottom, label + scale	clear, flat bottom	amber, flat bottom	clear, flat bottom, label + scale	amber, flat bottom, label + scale	clear, flat bottom	amber, flat bottom	
100 / PP box				100 / PP box				100 / PP box	100 / PP box	
To be used with*				To be used with*				To be used with*		
702068		702066		702032	702033	702034		702113	702051	702052
N 8 PP screw cap, black, closed top, Red Rubber / FEP colorless, 45° shore A, 1.3 mm		N 8 PP screw cap, black, closed top, Silicone white / PTFE red, 45° shore A, 1.3 mm		N 9 PP screw cap blue, closed top, PTFE virginal, white, 53° shore D, 0.25 mm	N 9 PP screw cap blue, closed top, Red Rubber / FEP colorless, 45° shore A, 1.0 mm	N 9 PP screw cap blue, closed top, Silicone white / PTFE red, 45° shore A, 1.0 mm		N 13 PP bonded screw cap, white, closed top, Silicone white / PTFE blue, 55° shore A, 1.3 mm	N 13 PP screw cap, black, closed top, Red Rubber / FEP colorless, 40° shore A, 1.5 mm	N 13 PP screw cap, black, closed top, Silicone white / PTFE red, 40° shore A, 1.3 mm
100 / PE bag		100 / PE bag		100 / PE bag	100 / PE bag	100 / PE bag		100 / PE bag	100 / PE bag	100 / PE bag

* Scale (1:1.6)

* Scale (1:1.6)

Container for screw neck vials N 15

REF	Description	Pack of
702518	36 position container blue with removable divider for screw neck vials N 15, 130 x 130 x 80 mm, with transparent lid (suitable for freezers)	1



■ Screw neck vials and screw closures for storage of liquid samples

N 15				N 18	N 20
15-425 thread				18-400 thread	20-400 thread
					
702096	702311	70285	702097	702098	702099
8 mL	8 mL	12 mL	12 mL	16 mL	24 mL
16.6 x 61 mm	16.6 x 61 mm	18.5 x 66 mm	18.5 x 66 mm	20.6 x 71 mm	22.7 x 86 mm
clear, flat bottom	amber, flat bot.	clear, flat bottom	amber, flat bot.	clear, flat bottom	clear, flat bottom
100 / PP box	100 / carton	100 / PP box	100 / carton	100 / PP box	100 / PP box
To be used with*				To be used with*	To be used with*
					
702114	702180			702115	702116 702181
N 15 PP bonded screw cap, white, closed top, Silicone white / PTFE blue, 55° shore A, 1.3 mm		N 15 PP bonded screw cap, black, center hole, Silicone white / PTFE beige, 45° shore A, 1.5 mm		N 18 PP bonded screw cap, white, closed top, Silicone white / PTFE blue, 55° shore A, 1.3 mm	N 20 PP bonded screw cap, white, closed top, Silicone white / PTFE blue, 55° shore A, 1.3 mm
					N 20 PP bonded screw cap, white, center hole, Silicone white / PTFE beige, 45° shore A, 1.5 mm
100 / PE bag	100 / PE bag			100 / PE bag	100 / PE bag 100 / PE bag

* Scale (1:1.6)

For screw neck vials with volumes larger than 24 mL please see page 43

DIV Diverse vials and caps

Snap cap vials (for storage of powdery samples) N 18 + N 22 with PE snap caps

N 18		N 22	
			
70271	70272	702019	70273
5 mL	10 mL	15 mL	25 mL
20 x 40 mm	22 x 50 mm	26 x 48 mm	26 x 65 mm
clear, 3 rd hydrol. class, flat bottom	clear, 3 rd hydrol. class, flat bottom	clear, 3 rd hydrol. class, flat bottom	clear, 3 rd hydrol. class, flat bottom
100 / carton	100 / carton	100 / carton	100 / carton
To be used with N 18 PE snap cap	To be used with N 18 PE snap cap	To be used with N 22 PE snap cap	To be used with N 22 PE snap cap
			
70274	70274	70275	70275
100 / PE bag	100 / PE bag	100 / PE bag	100 / PE bag



Safe handling of samples in headspace analysis

10 mL and 20 mL screw neck vials N 18 from MACHEREY-NAGEL with magnetic screw caps are designed for use on headspace autosamplers such as the CTC Combi PAL®.

They offer a more convenient, consistent and secure sealing process than crimping (better reproducibility). Due to the thinner septum (1.5 mm instead of 3 mm), needle penetration is easier and associated with less fragmentation (especially important for SPME applications). The screw caps provide a flat surface for the magnet; thus the filled vial cannot drop down from the magnet, which can be the case with overcrimped vials due to the convexity of the cap going along with it.

The closures are offered with a variety of different septa materials, with silicone/PTFE being recommended for its excellent analytical purity, high temperature resistance and outstanding penetration properties. Closed top closure variants are also available for sample storage / transport.



Headspace screw neck vials N 18

Precision thread, not standardized



702866	702866.2	702826	702826.0	702826.2
10 mL	10 mL	20 mL	20 mL	20 mL
22.5 x 46 mm	22.5 x 46 mm	22.5 x 75.5 mm	22.5 x 75.5 mm	22.5 x 75.5 mm
clear, rounded bottom	amber, rounded bottom	clear, rounded bottom	clear, with label, rounded bottom	amber, rounded bottom
100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box

Ready assembled, magnetic screw closures N 18

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702827		N 18 magnetic screw cap, silver, center hole	Silicone blue transparent / PTFE white	45° shore A	1.5 mm	100 / PE bag
702055		N 18 magnetic screw cap, silver, center hole	Silicone white / PTFE blue	45° shore A	1.5 mm	100 / PE bag
702140		as above, but with closed top				
702136		N 18 magnetic screw cap, silver, center hole	Silicone white / PTFE blue, slit	45° shore A	1.5 mm	100 / PE bag
702137		N 18 magnetic screw cap, silver, center hole	Butyl red / PTFE gray	50° shore A	1.5 mm	100 / PE bag
702139		as above, but with closed top				
702072		N 18 magnetic screw cap, silver, center hole	Red Rubber / TEF colorless	45° shore A	1.5 mm	100 / PE bag
702110		- (septa only)	Silicone white / PTFE blue	45° shore A	1.5 mm	100 / PE bag
702981		- (septa only)	Silicone blue transparent / PTFE white	45° shore A	1.5 mm	100 / PE bag
702138		- (septa only)	Butyl red / PTFE gray	50° shore A	1.5 mm	100 / PE bag

N 20 Crimp neck vials and caps N 20

Crimp neck vials N 20

5–10 mL



70204.36 70215.36	702917	702020	70205.36 70216.36	702918	702924
Varian®	PerkinElmer®	Metrohm®	Varian®	DANI, Agilent®	CTC, Varian®
5 mL	6 mL	5 mL	10 mL	10 mL	10 mL
20.5 x 38 mm	22 x 38.25 mm	21.7 x 38.25 mm	20.5 x 54.5 mm	22.5 x 46 mm	22.5 x 46 mm
clear / amber	clear	clear	clear / amber	clear	clear
flat bottom, flat DIN crimp neck	rounded bottom, bev. HS crimp neck	flat bottom, beveled HS crimp neck	flat bottom, flat DIN crimp neck	flat bottom, flat DIN crimp neck	rounded bottom, flat DIN crimp neck
100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box

20 mL



70206.36 70217.36	702086	70254	702540	702261	702263 702263.2	702263.0
		PerkinElmer®	PerkinElmer®	DANI, Agilent®	CTC, Varian®	CTC, Varian®
20 mL	20 mL	20 mL	20 mL	20 mL	20 mL	20 mL
23.25 x 75.5 mm	22 x 75 mm	23 x 75.5 mm	23 x 75.5 mm	22.5 x 75.5 mm	22.5 x 75.5 mm	22.5 x 75.5 mm
clear / amber	clear, 15 mL graduation	clear	clear, with label	clear	clear / amber	clear, with label
flat bottom, flat DIN crimp neck	rounded bottom, bev. HS crimp neck	rounded bottom, bev. HS crimp neck	rounded bottom, bev. HS crimp neck	flat bottom, flat DIN crimp neck	rounded bottom, flat DIN crimp neck	rounded bottom, flat DIN crimp neck
100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box	100 / PP box

Crimp neck vials N 20

25 mL, 50 mL and 100 mL



70210.36	70208.36	70209.1
25 mL	50 mL	100 mL
30 x 65 mm	31 x 101 mm	51.6 x 94.5 mm
clear	clear	clear, 3 rd hydrolytic class
flat bottom	flat bottom	flat bottom
flat DIN crimp neck	flat DIN crimp neck	flat DIN crimp neck
100 / PP box	100 / PP box	88 / shrink-wrapped tray

Crimping tools N 20

REF	Type of crimping tool	Pack of
735120	Manual crimper (standard), height adjustable, for 20 mm aluminum crimp caps	1 / box
735132	Manual crimper (standard), height adjustable, for 20 mm Flip Top/Flip Off caps	1 / box
735920	Manual decapper (standard) for 20 mm aluminum crimp caps	1 / box
735220	Ergonomic crimper for 20 mm aluminum crimp caps	1 / box
735320	Ergonomic decapper for 20 mm aluminum crimp caps	1 / box
735520	Electronic crimper for 20 mm aluminum crimp caps (battery-powered)	1 / box
735532	Electronic crimper for 20 mm Flip Top/Flip Off caps (battery-powered)	1 / box
735620	Electronic decapper for 20 mm aluminum crimp caps (battery-powered)	1 / box
735700	Electronic high power crimping tool with power supply	1 / box
735720	Crimping head for 20 mm crimp caps (aluminum, magnetic, bi-metal)	1 / box
735732	Crimping head for 20 mm Flip Top/Flip Off caps	1 / box
735820	Decapping head for 20 mm crimp caps (aluminum, magnetic, bi-metal)	1 / box



N 20 Crimp neck vials and caps N 20

Ready assembled crimp closures, plain crimp caps and single septa N 20

REF	(Scale 1:2)	Cap description	Septum description	Hardness	Thickness	Pack of
Center hole caps						
702773		N 20 aluminum crimp cap, silver, center hole	Butyl red / PTFE gray	50° shore A	3 mm	100 / PE bag
702775		N 20 aluminum crimp cap, silver, center hole	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100 / PE bag
70234.9		N 20 aluminum crimp cap, silver, center hole	Molded septum Butyl / PTFE gray	50° shore A	3 mm	100 / PE bag
70234		N 20 aluminum crimp cap, silver, center hole	Butyl dark gray / PTFE gray*	50° shore A	3 mm	100 / PE bag
70234.10		N 20 aluminum crimp cap, silver, center hole	Butyl dark gray / PTFE gray* (high purity)	50° shore A	3 mm	100 / PE bag
70237.2		N 20 aluminum crimp cap, silver, center hole	Bromobutyl stopper, gray, unsembled (separate parts)	50° shore A		100 each / PE bags
702093		N 20 aluminum crimp cap, silver, center hole	Silicone blue transparent / PTFE colorless	40° shore A	3 mm	100 / PE bag
702340		N 20 aluminum crimp cap, silver, center hole	Silicone blue transparent / central PTFE lamination colorless*	45° shore A	3 mm	100 / PE bag
702094		N 20 aluminum crimp cap, silver, center hole	Silicone white / PTFE beige	40° shore A	3 mm	100 / PE bag
702145		N 20 aluminum crimp cap, silver, center hole	Silicone white / FEP-/Aluminum foil silver	45° shore A	3 mm	100 / PE bag
702804		N 20 aluminum crimp cap, silver, center hole (no liner)	-	-	-	100 / PE bag
702112		N 20 aluminum crimp cap, gold, center hole (no liner)	-	-	-	100 / PE bag
Pressure release caps						
702829		N 20 aluminum pressure release cap, silver, center hole	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100 / PE bag
70234.8		N 20 aluminum pressure release cap, silver, center hole	Molded septum Butyl / PTFE gray	50° shore A	3 mm	100 / PE bag
702071		N 20 aluminum pressure release cap, silver, center hole	Butyl dark gray / PTFE gray* (high purity)	50° shore A	3 mm	100 / PE bag
702927		N 20 aluminum pressure release cap, silver, center hole	Silicone blue transparent / PTFE colorless	45° shore A	3 mm	100 / PE bag
702835		N 20 aluminum pressure release cap, silver, center hole	Silicone white / PTFE beige	45° shore A	3 mm	100 / PE bag
702799		N 20 aluminum pressure release cap, silver, center hole (no liner)	-	-	-	100 / PE bag

* only centrally laminated with PTFE, typically called Pharma-Fix

 Ready assembled crimp closures, plain crimp caps and single septa N 20

REF	(Scale 1:2)	Cap description	Septum description	Hardness	Thickness	Pack of
 Bi-metal crimp caps						
702838		N 20 bi-metal crimp cap, blue/silver, center hole	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100 / PE bag
702834		N 20 bi-metal crimp cap, blue/silver, center hole	Silicone blue transparent / PTFE colorless	40° shore A	3 mm	100 / PE bag
702341		N 20 bi-metal crimp cap, blue/silver, center hole	Silicone blue transparent / central PTFE lamination colorless*	45° shore A	3 mm	100 / PE bag
702837		N 20 bi-metal crimp cap, blue/silver, center hole	Silicone white / PTFE beige	40° shore A	3 mm	100 / PE bag
702833		N 20 bi-metal crimp cap, blue/silver, center hole (no liner)	-	-	-	100 / PE bag
 Magnetic crimp caps						
702928		N 20 magnetic crimp cap, silver, 8 mm center hole	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100 / PE bag
702928.9		N 20 magnetic crimp cap, silver, 8 mm center hole	Butyl dark gray / PTFE gray*	50° shore A	3 mm	100 / PE bag
702929		N 20 magnetic crimp cap, silver, 8 mm center hole	Silicone blue transparent / PTFE colorless	40° shore A	3 mm	100 / PE bag
702808		N 20 magnetic crimp cap, silver, 8 mm center hole (no liner)	-	-	-	100 / PE bag
 Center tear off caps						
70233		N 20 aluminum center tear off cap, gold	 Butyl dark gray / PTFE gray*	50° shore A	3 mm	100 / PE bag
70236.1		N 20 aluminum center tear off cap, silver (no liner)	-	-	-	100 / PE bag
 Complete tear off caps						
70235		N 20 aluminum complete tear off cap, silver	Butyl dark gray / PTFE gray*	50° shore A	3 mm	100 / PE bag
702839		N 20 aluminum complete tear off cap, silver	Silicone white / PTFE beige	40° shore A	3 mm	100 / PE bag
702805		N 20 aluminum complete tear off cap, silver (no liner)	-	-	-	100 / PE bag
 Single septa						
70277		- (septum only)	Butyl red / PTFE gray	50° shore A	3 mm	100 / PE bag
702057		- (septum only)	Butyl light gray / PTFE dark gray	50° shore A	3 mm	100 / PE bag
702101		- (septum only)	Molded septum Butyl / PTFE gray	50° shore A	3 mm	100 / PE bag
702D20TB		- (septum only)	Butyl dark gray / PTFE gray*	50° shore A	3 mm	100 / PE bag
702780		- (septum only)	Silicone blue transparent / PTFE colorless	40° shore A	3 mm	100 / PE bag
70278		- (septum only)	Silicone white / PTFE beige	40° shore A	3 mm	100 / PE bag

* only centrally laminated with PTFE, typically called Pharma-Fix

N 20 Crimp neck vials and caps N 20

Stoppers N 20

REF	(Scale 1:2)	Cap description	Description	Hardness	Thickness	Pack of
702931.3		- (stopper only)	N 20 stopper, Bromobutyl gray	50° shore A	–	100 / PE bag
702931.1		- (stopper only)	N 20 stopper, Bromobutyl red	45° shore A	–	100 / PE bag

PE caps N 20 and suitable single septa N 19

REF	(Scale 1:2)	Cap description	Septum description	Hardness	Thickness	Pack of
Height: 8.4 mm → for vials with beveled HS crimp neck						
70266		N 20 PE cap, transparent, for beveled HS crimp neck N 20, 4.3 mm center hole	(no liner)	-	-	100 / PE bag
70242		cap as above	Butyl beige / PTFE gray unassembled (as separate parts)	55° shore A	1.3 mm	100 each / PE bags
702128		cap as above	Natural rubber red-orange / TEF colorless, ready assembled	45° shore A	1.3 mm	100 / PE bag
Height: 9.1 mm → for vials with flat DIN crimp neck						
70267		N 20 PE cap, transparent, for flat DIN crimp neck N 20, 4.3 mm center hole	(no liner)	-	-	100 / PE bag
70240		cap as above	Butyl beige / PTFE gray unassembled (as separate parts)	55° shore A	1.3 mm	100 each / PE bags
702129		cap as above	Natural rubber red-orange / TEF colorless, ready assembled	45° shore A	1.3 mm	100 / PE bag
70269		- (septum only)	Butyl beige / PTFE gray	55° shore A	1.3 mm	100 / PE bag
702904		- (septum only)	Natural rubber red-orange / TEF colorless	45° shore A	1.3 mm	100 / PE bag
702144		- (septum only)	Silicone blue transparent / PTFE white	45° shore A	1.3 mm	100 / PE bag

Container for screw neck vials N 18 and crimp neck vials N 20

REF	Description	Pack of
702516	25 position container blue for 10 mL + 20 mL screw neck vials N 18 and crimp neck vials N 20, 130 x 130 x 80 mm, with transparent lid (suitable for freezers)	1



Vial rack for screw neck vials N 18 and crimp neck vials N 20

REF	Description	Pack of
702503	36 position polypropylene vial rack blue, coded, for vials with a diameter of 23.1 mm max. Dimensions: 323 x 91 x 30 mm, stackable	5 / carton



Screw neck vials N 24 (EPA)

For TOC / VOC analysis and sample storage, 24-414 thread

			
702021 / 702022	702132 / 702133	702023 / 702024	702074 / 702131
20 mL	30 mL	40 mL	60 mL
27.5 x 57 mm	27.5 x 72.5 mm	27.5 x 95 mm	27.5 x 140 mm
clear / amber, flat bottom	clear / amber, flat bottom	clear / amber, flat bottom	clear / amber, flat bottom
100 / PP box	100 / PP box	100 / PP box	100 / PP box







Container for screw neck vials N 24

REF	Description	Pack of
702517	16 position container blue, with removable divider for screw neck vials N 24 (20 mL, 30 mL, 40 mL) 130 x 130 x 102 mm, with transparent lid (suitable for freezers)	1



N 24 Screw neck vials and caps N 24

Bonded and ready assembled screw closures + plain screw caps N 24, single septa N 22

REF	(Scale 1:1.6)	Cap description	Septum description	Hardness	Thickness	Pack of
702058		N 24 PP bonded screw cap, white, center hole	Silicone white / PTFE beige	45° shore A	3.2 mm	100 / PE bag
702059		N 24 PP bonded screw cap, white, closed top	Silicone white / PTFE beige	45° shore A	3.2 mm	100 / PE bag
702073		N 24 PP bonded screw cap, white, center hole	Red Rubber / TEF colorless	45° shore A	2.5 mm	100 / PE bag
702168		N 24 PP bonded screw cap, white, center hole, with removable protection lid	Silicone natural / PTFE colorless	45° shore A	3.2 mm	100 / PE bag
702130		N 24 PP screw cap, white, center hole	Butyl red / PTFE gray	50° shore A	2.4 mm	100 / PE bag
702102		N 24 PP screw cap, white, closed top	Butyl red / PTFE gray	50° shore A	2.4 mm	100 / PE bag
702171		N 24 PP screw cap, white, center hole	Silicone white / PTFE beige	45° shore A	3.2 mm	100 / PE bag
702060		N 24 PP screw cap, white, center hole	no liner	–	–	100 / PE bag
702061		N 24 PP screw cap, white, closed top	no liner	–	–	100 / PE bag
702062		- (septa only)	Silicone natural / PTFE colorless	45° shore A	3.2 mm	100 / PE bag
702791		- (septa only)	Butyl red / PTFE gray	50° shore A	2.4 mm	100 / PE bag

Pre-sealed vial-closure combinations with screw neck N 24

Vial description	Closure description	Pack of	REF
Pre-sealed vials 702021 : 20 mL screw neck vial N 24, 27.5 x 57 mm, clear, flat bottom	pre-screwed with 702073 : N 24 PP screw cap (bonded), white, center hole, Red Rubber / TEF colorless, 45° shore A, 2.5 mm	100	702865
Pre-sealed vials 702021 : 20 mL screw neck vial N 24, 27.5 x 57 mm, clear, flat bottom	pre-screwed with 702058 : N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702894
Pre-sealed vials 702021 : 20 mL screw neck vial N 24, 27.5 x 57 mm, clear, flat bottom	pre-screwed with 702059 : N 24 PP screw cap (bonded), white, closed top, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702884
Pre-sealed vials 702132 : 30 mL screw neck vial N 24, 27.5 x 72.5 mm, clear, flat bottom	pre-screwed with 702058 : N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702831
Pre-sealed vials 702023 : 40 mL screw neck vial N 24, 27.5 x 95 mm, clear, flat bottom	pre-screwed with 702058 : N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702877
Pre-sealed vials 702023 : 40 mL screw neck vial N 24, 27.5 x 95 mm, clear, flat bottom	pre-screwed with 702168 : N 24 PP screw cap (bonded), white, 15 mm center hole, with removable protection lid, Silicone natural / PTFE colorless, 45° shore A, 3.2 mm	100	702872
Pre-sealed vials 702131 : 60 mL screw neck vial N 24, 27.5 x 140 mm, amber, flat bottom	pre-screwed with 702058 : N 24 PP screw cap (bonded), white, center hole, Silicone white / PTFE beige, 45° shore A, 3.2 mm	100	702886
Other pre-sealed vial-closure combinations on request			



Universally usable screw closure for storage and analysis

Explore our video on this unique and innovative closure (MN REF 702168) with removable protection lid for the underlying center hole that allows injection and analysis after storage without the need to change the cap. The lid keeps the injection area clean from dust or any other contamination. The sealing septa consist of high purity Silicone/PTFE being bonded to the cap ensuring safe penetration with a needle.

Watch the video under <https://youtu.be/vwHnb-k8QJI>



Vial accessories and other consumables for chromatography

GC injection port septa

Premium HT-quality – Blister packs for cleanliness and easy handling



Septum type		BTO septa	AG3 septa	Marathon septa
				
OD	T _{max}	for highest demands in GC and GC-MS – CenterGuide™ –	with higher durability than BTO – CenterGuide™ –	with extreme durability > 400 injections – CenterGuide™ –
9 mm	400 °C	702646	702656	702660
11 mm	400 °C	702647	702657	702661
11.5 mm	400 °C	702648	702658	702662
Shimadzu®	300 °C	702649	702659	702663
Pack of		25 / blister	25 / blister	25 / blister

Standard quality

Septum type	Standard septa (ST)	High temperature septa (HT)	Silicone septa	Silicone / PTFE septa
OD	beige silicone, 60° shore A, 4 mm	red, non-bleeding silicone, 60° shore A, 3 mm (320 °C max.)	soft, transparent 3 mm (250 °C max.)	white silicone, PTFE 3 mm (200 °C max.)
9 mm	702609	702619	702602	
10 mm	702610	702620		702625
11 mm	702611	702621	702604	702626
12 mm	702612	702622	702605	702627
13 mm	702613	702623	702606	702628
17 mm		702632		
Pack of	50 / plastic container	50 / plastic container	50 / plastic container	50 / plastic container

Capillary connectors for GC

REF	Capillary connector		Pack of
707971	Universal, linear capillary connector for 2 columns with an ID of 0.20–0.53 mm		5 / box
707972	Universal, linear capillary connector for 2 columns with an ID of 0.20–0.53 mm		10 / box
707973	Universal, linear capillary connector for 3 columns with an ID of 0.20–0.53 mm, Y splitter		1 / box

Manual crimping tools – Standard

User profile

Manual crimping tool - Standard





Cap sizes (N 8, N 11, N 13, N 20)

Aluminum

Bi-metal

Magnetic*



Number of crimp processes

Reproducibility

Durability

Mobility

Comfort adjust.

Convenience



REF	Description	Pack of
Crimpers		
735126	Manual crimper (standard), for 8 mm crimp caps	1 / box
735111	Manual crimper (standard), height adjustable, for 11 mm crimp caps	1 / box
735113	Manual crimper (standard), height adjustable, for 13 mm crimp caps	1 / box
735133	Manual crimper (standard), height adjustable, for 13 mm Flip Top/Flip Off caps	1 / box
735120	Manual crimper (standard), height adjustable, for 20 mm crimp caps	1 / box
735132	Manual crimper (standard), height adjustable, for 20 mm Flip Top/Flip Off caps	1 / box
Decappers		
735408	Manual decapper (standard), for 8 mm crimp caps	1 / box
735911	Manual decapper (standard), for 11 mm crimp caps	1 / box
735913	Manual decapper (standard), for 13 mm crimp caps	1 / box
735920	Manual decapper (standard), for 20 mm crimp caps	1 / box



Vial accessories and other consumables for chromatography

Manual crimping tools – Ergonomic

User profile


Manual crimping tool - Ergonomic



Cap sizes (N 8, N 11, N 20)

Aluminum

Bi-metal



Required force

Number of crimp processes

Reproducibility

Durability

Mobility

Comfort adjustments

Convenience



REF	Description	Pack of
Crimpers		
735208	Manual crimper (ergonomic), adjustable, for 8 mm crimp caps	1 / box
735211	Manual crimper (ergonomic), adjustable, for 11 mm crimp caps	1 / box
735220	Manual crimper (ergonomic), adjustable, for 20 mm crimp caps	1 / box
Decappers		
735311	Manual decapper (ergonomic), for 11 mm crimp caps	1 / box
735320	Manual decapper (ergonomic), for 20 mm crimp caps	1 / box
Accessories for ergonomic crimping tools		
735509	Rack for two ergonomic crimping tools, blue, 240 x 95 x 65 mm	1

Electronic crimping tools – Battery-powered

User profile

Battery-powered – electr. crimping tool



Cap sizes (N 11, N 13, N 20)

Aluminum



Required force

Number of crimp processes

Reproducibility

Durability

Mobility

Comfort adjustments

Convenience




REF	Description	Pack of
Crimpers		
735511	Electronic crimper (battery-powered), for 11 mm aluminum crimp caps	1 / box
735513	Electronic crimper (battery-powered), for 13 mm aluminum crimp caps	1 / box
735533	Electronic crimper (battery-powered), for 13 mm Flip Top/Flip Off caps	1 / box
735520	Electronic crimper (battery-powered), for 20 mm aluminum crimp caps	1 / box
735532	Electronic crimper (battery-powered), for 20 mm Flip Top/Flip Off caps	1 / box
Decappers		
735611	Electronic decapper (battery-powered), for 11 mm aluminum crimp caps	1 / box
735613	Electronic decapper (battery-powered), for 13 mm aluminum crimp caps	1 / box
735620	Electronic decapper (battery-powered), for 20 mm aluminum crimp caps	1 / box
Accessories for battery-powered, electronic crimping tools		
735500	Battery 6.6 Volt, 8.6 Wh	1 / box
735509	Rack for two electronic crimping tools, blue, 240 x 95 x 65 mm	1

Vial accessories and other consumables for chromatography

Electronic crimping tools – High Power crimping tool

User profile

Electronic high power crimping tool



Required force

Cap sizes (N 8, N 11, N 13, N 20)

Aluminum

Bi-metal

Magnetic

Number of crimp processes

Reproducibility

Durability

Mobility*

Comfort adjustments

Convenience



* normally stationary usage with a stand

REF	Description	Pack of
Basic high power crimping tool		
735700	Electronic high power crimping tool with fixed power supply	1 / box
Crimping heads		
735708	Crimping head for 735700 for 8 mm crimp caps	1 / box
735711	Crimping head for 735700 for 11 mm crimp caps	1 / box
735713	Crimping head for 735700 for 13 mm crimp caps	1 / box
735733	Crimping head for 735700 for 13 mm Flip Top/Flip Off caps	1 / box
735720	Crimping head for 735700 for 20 mm crimp caps	1 / box
735732	Crimping head for 735700 for 20 mm Flip Top/Flip Off caps	1 / box
Decapping heads		
735811	Decapping head for 735700 for 11 mm crimp caps	1 / box
735813	Decapping head for 735700 for 13 mm crimp caps	1 / box
735820	Decapping head for 735700 for 20 mm crimp caps	1 / box
Accessories for high power crimping tool		
735501	Stand for electronic crimping tools	1 / box
735509	Rack for two electronic HP crimping tools, blue, 240 x 95 x 65 mm	1

Accessories for crimping tools

REF	Description	Pack of
735500	Battery 6.6 Volt, 8.6 Wh	1
735501	Stand for electronic crimping tools	1
735509	Rack for two crimping tools either type ergonomic (manual) or electronic, blue, 240 x 95 x 65 mm	1



Vial containers

REF	Description	Pack of
702514	81 position container blue, with firmly integrated divider, for all vials 11.6 x 32 mm, 130 x 130 x 45 mm, coded, with transparent lid (suitable for freezers)	1
702515	49 position container blue, with firmly integrated divider, for all vials N 13 with flat bottom, 130 x 130 x 50 mm, with transparent lid (suitable for freezers)	1
702516	25 position container blue, with removable divider, for all 10 mL and 20 mL screw neck vials N 18 and crimp neck vials N 20, 130 x 130 x 80 mm, with transparent lid (suitable for freezers)	1
702518	36 position container blue, with removable divider, for screw neck vials N 15 (sample storage: 702096, 702311, 70285, 702097) 130 x 130 x 80 mm, with transparent lid (suitable for freezers)	1
702517	16 position container blue, with removable divider, for all 20 mL, 30 mL and 40 mL screw neck vials N 24, 130 x 130 x 102 mm, with transparent lid (suitable for freezers)	1



Vial racks

REF	Description	Pack of
702502	50 position polypropylene vial rack blue, for all vials 11.6 x 32 mm with flat bottom Dimensions: 190 x 100 x 22 mm, stackable	1
702504	50 position polypropylene vial rack blue, for all vials with a diameter of 15 mm max. and flat bottom Dimensions: 240 x 120 x 28 mm, stackable	1
702503	36 position polypropylene vial rack blue, for all vials with a diameter of 23.1 mm max. and flat bottom Dimensions: 323 x 91 x 30 mm, stackable	5 / carton



Numerical index

REF	Product	Page
702D20TB	 septum N 20	41
702001	 crimp closure N 11	28
702004	 screw neck vial N 8	19
	 screw neck vial N 8	34
702005	 insert for vials N 8	19
	 insert for vials N 11, small opening	26
702006	 screw neck vial N 9	21
702007	 screw neck vial N 9	21
702007.0	 screw neck vial N 9	21
702008	 screw neck vial N 9	21
702009	 screw neck vial N 9	21
70201CG	 crimp neck vial N 11	26
70201HP	 crimp neck vial N 11	26
70201HP.2	 crimp neck vial N 11	26
702010	 screw neck vial N 9	21
702011	 screw neck vial N 10	25
702012	 screw neck vial N 10	25
702013	 screw neck vial N 10	25
702014	 crimp neck vial N 11	27
702015	 crimp neck vial N 11	27
702016	 crimp neck vial N 11	27
702019	 snap cap vial N 22	36
70202.1	 shell vial N 8	18
702020	 crimp neck vial N 20	38
702021	 screw neck vial N 24	43
702022	 screw neck vial N 24	43
702023	 screw neck vial N 24	43
702024	 screw neck vial N 24	43
702025	 crimp closure N 8	18
702026	 bonded screw closure N 9	23
702027	 bonded screw closure N 9	23
702028	 bonded screw closure N 9	23
702029	 screw closure N 9	22
70203	 crimp neck vial N 13	32
702030	 screw closure N 9	22
702031	 screw closure N 9	22
702032	 screw closure N 9	22
	 screw closure N 9	34
702033	 screw closure N 9	22
	 screw closure N 9	34
702034	 screw closure N 9	22
	 screw closure N 9	34
702035	 screw closure N 9	22
702036	 screw closure N 9	22
702037	 screw closure N 9	22
702038	 screw closure N 9	22
702039	 screw closure N 9	22
70204.36	 crimp neck vial N 20	38
702040	 screw closure N 9	22
702041	 septum N 9	23
702042	septum N 9	23
702043	septum N 9	23
702044	bonded screw closure N 10	25
702046	bonded screw closure N 10	25

REF	Product	Page
702047	 bonded screw closure N 10	25
702048	 screw closure N 10	25
702049	 screw cap N 10	25
70205.36	 crimp neck vial N 20	38
702050	 screw closure N 13	33
702051	 screw closure N 13	33
	 screw closure N 13	34
702052	 screw closure N 13	33
	 screw closure N 13	34
702053	 septum N 12	33
702055	 magnetic screw closure N 18	37
702057	 septum N 20	41
702058	 bonded screw closure N 24	44
702059	 bonded screw closure N 24	44
70206.36	 crimp neck vial N 20	38
702060	 screw cap N 24	44
702061	 screw cap N 24	44
702062	 septum N 22	44
702063	 snap ring closure N 11	31
702063.2080	 snap ring closure N 11	31
702064	 snap ring closure N 11	31
702066	 screw closure N 8	20
	 screw closure N 8	34
702067	 screw closure N 8	20
702068	 screw closure N 8	20
	 screw closure N 8	34
702069	 screw closure N 8	20
702070	 septum N 8	20
702071	 crimp closure N 20	40
702072	 magnetic screw closure N 18	37
702073	 bonded screw closure N 24	44
702074	 screw neck vial N 24	43
702075	 crimp neck vial N 11	26
702076	 crimp neck vial N 11	26
702077	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
702078	 screw neck vial N 9	21
	 screw neck vial N 9	34
702079	 screw neck vial N 9	21
	 screw neck vial N 9	34
70208.36	 crimp neck vial N 20	39
702080	 screw closure N 9	22
702081	 screw closure N 9	22
702082	 screw closure N 9	22
702083	 screw closure N 9	22
702084	 screw closure N 9	22
702085	 screw closure N 9	22
702086	 crimp neck vial N 20	38
702088	 screw neck vial N 9	21
70209.1	crimp neck vial N 20	39
702093	crimp closure N 20	40
702094	crimp closure N 20	40
702095	snap ring closure N 11	31

Numerical index

REF	Product	Page
702096	 screw neck vial N 15	35
702097	 screw neck vial N 15	35
702098	 screw neck vial N 18	35
702099	 screw neck vial N 20	35
70210.36	 crimp neck vial N 20	39
702101	 septum N 20	41
702101HP	 pre-sealed vial-closure combination N 11	29
702102	 screw closure N 24	44
702107	 screw closure N 9	22
702108	 snap ring closure N 11	31
702109	 screw closure N 9	22
702110	 septum N 18	37
702112	 crimp cap N 20	40
702113	 bonded screw closure N 13	33
	 bonded screw closure N 13	34
702114	 bonded screw closure N 15	35
702115	 bonded screw closure N 18	35
702116	 bonded screw closure N 20	35
702128	 PE cap N 20	42
702129	 PE cap N 20	42
70213	 screw neck vial N 8	19
	 screw neck vial N 8	34
70213.2	 screw neck vial N 8	19
	 screw neck vial N 8	34
702130	 screw closure N 24	44
702131	 screw neck vial N 24	43
702132	 screw neck vial N 24	43
702133	 screw neck vial N 24	43
702134	 crimp neck vial N 11	27
	 snap ring vial N 11	30
702135	 screw neck vial N 9	21
702136	 magnetic screw closure N 18	37
702137	 magnetic screw closure N 18	37
702138	 septum N 18	37
702139	 magnetic screw closure N 18	37
70214CG	 crimp neck vial N 11	26
702140	 magnetic screw closure N 18	37
702141	 crimp neck vial N 11	27
702142	 snap ring closure N 11	31
702143	 snap ring closure N 11	31
702144	 septum N 19	42
702145	 crimp closure N 20	40
702147	 screw closure N 9	22
702148	 septum N 9	23
702150	 snap ring closure N 11	31
702151	 snap ring closure N 11	31
70215.36	 crimp neck vial N 20	38
702155	 magnetic screw closure N 9	23
702156	 magnetic screw closure N 9	23
702158	 screw closure N 9	22
702159	 screw closure N 9	22
70216.36	 crimp neck vial N 20	38
702160	 screw cap N 9	22
702161	screw cap N 9	22

REF	Product	Page
702162	 screw cap N 9	22
702163	 screw cap N 9	22
702164	 screw cap N 9	22
702165	 screw cap N 9	22
702168	 bonded screw closure N 24	44
70217.36	 crimp neck vial N 20	38
702170	 vial-insert combination N 11	30
702171	 screw closure N 24	44
702172	 screw neck vial N 9	21
702173	 snap ring vial N 11	30
702174	 snap ring vial N 11	30
702175	 crimp closure N 11	28
702176	 vial-insert combination N 11	30
702177	 vial-insert combination N 9	21
702178	 vial-insert combination N 9	21
702179	 vial-insert combination N 9	21
702180	 bonded screw closure N 15	35
702181	 bonded screw closure N 20	35
702261	 crimp neck vial N 20	38
702263	 crimp neck vial N 20	38
702263.0	 crimp neck vial N 20	38
702263.2	 crimp neck vial N 20	38
702282	 screw neck vial N 9	21
	 screw neck vial N 9	34
702283	 screw neck vial N 9	21
	 screw neck vial N 9	34
702284	 screw neck vial N 9	21
	 screw neck vial N 9	34
702286	 screw closure N 9	22
702287	 screw closure N 9	22
702287.1	 screw closure N 9	22
702288	 screw closure N 9	22
702288.1	 screw closure N 9	22
702292	 septum N 12	33
702293	 screw neck vial N 9	21
	 screw neck vial N 9	34
702295	 snap ring closure N 11	31
70231	 crimp closure N 11	28
70231.1	 crimp closure N 11	28
70231.2	 crimp closure N 11	28
702311	 screw neck vial N 15	35
70232	 crimp closure N 13	32
70233	 crimp closure N 20	41
702334	 crimp neck vial N 11	27
	 snap ring vial N 11	30
702335	 screw neck vial N 9	21
70234	 crimp closure N 20	40
70234.8	 crimp closure N 20	40
70234.9	 crimp closure N 20	40
70234.10	 crimp closure N 20 (high purity)	40
702340	 crimp closure N 20	40
702341	 crimp closure N 20	41
70235	 crimp closure N 20	41
70236.1	crimp cap N 20	41

Numerical index

REF	Product	Page
70237.2	 crimp closure N 20	40
70239	 crimp closure N 11	28
70240	 PE cap N 20	42
702401	 PE cap N 11	29
	 PE cap N 11	31
702402	 screw closure N 9	22
702403	 snap ring closure N 11	31
702405	 screw closure N 9	22
702406	 screw closure N 9	22
70242	 PE cap N 20	42
702437	 screw closure N 8	20
70245	 screw closure N 8	20
70246	 crimp closure N 11	28
70246.1	 crimp closure N 11	28
70246.2	 crimp closure N 11	28
70246.3	 crimp closure N 11	28
70248	 septum N 8	20
702481	 septum N 8	20
70249	 screw cap N 8	20
70250	 screw cap N 8	20
702500	 screw neck vial N 9	21
702502	 vial rack for vials N 8	19
	 vial rack for vials N 9	22
	 vial rack for vials N 10	25
	 vial rack for vials N 11	26
	 vial rack for vials N 8, N 9, N 10, N 11	50
702503	 vial rack for vials N 18, N 20	42
	 vial rack for vials N 18, N 20	50
702504	 vial rack for vials N 13	32
	 vial rack for vials N 13	50
702514	 container for vials N 8	19
	 container for vials N 9	22
	 container for vials N 10	25
	 container for vials N 11	26
	 container for vials N 8, N 9, N 10 and N 11	50
702515	 container for vials N 13	32
	 container for vials N 13	50
702516	 container for vials N 18, N 20	42
	 container for vials N 18, N 20	50
702517	 container for vials N 24	43
	 container for vials N 24	50
702518	 container for vials N 15	34
	 container for vials N 15	50
70254	 crimp neck vial N 20	38
702540	 crimp neck vial N 20	38
70256	 crimp closure N 11	28
70257	 crimp closure N 13	32
70260	 septum N 12	33
702602	 GC injection port septum	46
702604	 GC injection port septum	46
702605	GC injection port septum	46
702606	GC injection port septum	46
702609	GC injection port septum	46
70261	septum N 8	20

REF	Product	Page
702610	 GC injection port septum	46
702611	 GC injection port septum	46
702612	 GC injection port septum	46
702613	 GC injection port septum	46
702619	 GC injection port septum	46
702620	 GC injection port septum	46
702621	 GC injection port septum	46
702622	 GC injection port septum	46
702623	 GC injection port septum	46
702625	 GC injection port septum	46
702626	 GC injection port septum	46
702627	 GC injection port septum	46
702628	 GC injection port septum	46
702632	 GC injection port septum	46
702646	 GC injection port septum	46
702647	 GC injection port septum	46
702648	 GC injection port septum	46
702649	 GC injection port septum	46
702656	 GC injection port septum	46
702657	 GC injection port septum	46
702658	 GC injection port septum	46
702659	 GC injection port septum	46
70266	 PE cap N 20	42
702660	 GC injection port septum	46
702661	 GC injection port septum	46
702662	 GC injection port septum	46
702663	 GC injection port septum	46
70267	 PE cap N 20	42
70269	 septum N 19	42
702708	 snap ring vial N 11	30
702709	 snap ring vial N 11	30
70271	 snap cap vial N 18	36
702710	 snap ring closure N 11	31
702710.1	 snap ring closure N 11	31
702710.2080	 snap ring closure N 11	31
702712	 snap ring vial N 11	30
702713	 snap ring vial N 11	30
702714	 snap ring vial N 11	30
702716	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
702717.2	 snap ring closure N 11	31
702717.2080	 snap ring closure N 11	31
702718	 snap ring closure N 11	31
702718.1	 snap ring closure N 11	31
702718.2080	 snap ring closure N 11	31
70272	 snap cap vial N 18	36
70273	 snap cap vial N 22	36
702730	 crimp closure N 11	28
702730.1	 crimp closure N 11	28
702730.2	crimp closure N 11	28
702730.3	crimp closure N 11	28
702731	snap ring closure N 11	31
702732	screw closure N 9	22

Numerical index

REF	Product	Page
702736	 crimp closure N 11	28
70274	 PE snap cap N 18	36
70275	 PE snap cap N 22	36
70277	 septum N 20	41
702773	 crimp closure N 20	40
702775	 crimp closure N 20	40
70278	 septum N 20	41
702780	 septum N 20	41
702791	 septum N 22	44
702799	 crimp cap N 20	40
702800	 crimp cap N 8	18
702801	 crimp cap N 11	29
702802	 crimp cap N 13	32
702803	 crimp cap N 13	32
702804	 crimp cap N 20	40
702805	 crimp cap N 20	41
702807	 PE plug N 8	18
702808	 crimp cap N 20	41
702809	 snap ring vial N 11	30
702813	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
702818	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
702818.1	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
702819	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
70282	 crimp neck vial N 8	18
702823	 crimp closure N 11	28
702824	 insert for vials N 8	19
	 insert for vials N 11, small opening	26
702825	 insert for vials N 9	21
	 insert for vials N 10	25
	 insert for vials N 11, wide opening	26, 30
702826	 Headspace screw neck vial N 18	37
702826.0	 Headspace screw neck vial N 18	37
702826.2	 Headspace screw neck vial N 18	37
702827	 magnetic screw closure N 18	37
702829	 crimp closure N 20	40
702831	 pre-sealed vial-closure combination N 24	45
702833	 crimp cap N 20	41
702834	 crimp closure N 20	41
702835	 crimp closure N 20	40
702837	 crimp closure N 20	41
702838	 crimp closure N 20	41
702839	 crimp closure N 20	41
70284	 crimp closure N 11	28
70285	 screw neck vial N 15	35
702857	 pre-sealed vial-closure combination N 9	24
702858	pre-sealed vial-closure combination N 9	24

REF	Product	Page
702859	 pre-sealed vial-closure combination N 11	29
70286	 crimp neck vial N 8	18
702860	 screw neck vial N 8	19
702863	 pre-sealed vial-closure combination N 9	24
702864	 pre-sealed vial-closure combination N 9	24
702865	 pre-sealed vial-closure combination N 24	45
702866	 Headspace screw neck vial N 18	37
702866.2	 Headspace screw neck vial N 18	37
702867	 pre-sealed vial-closure combination N 11	29
702872	 pre-sealed vial-closure combination N 24	45
702873	 pre-sealed vial-closure combination N 9	24
702874	 pre-sealed vial-closure combination N 9	24
702877	 pre-sealed vial-closure combination N 24	45
702878	 crimp closure N 8	18
702879	 crimp closure N 11	28
70288	 crimp closure N 11	28
70288.1	 crimp closure N 11	28
70288.2	 crimp closure N 11	28
70288.3	 crimp closure N 11	28
702882	 pre-sealed vial-closure combination N 11	29
702884	 pre-sealed vial-closure combination N 24	45
702885	 crimp neck vial N 11	26
702886	 pre-sealed vial-closure combination N 24	45
702887	 pre-sealed vial-closure combination N 11	29
702888	 crimp neck vial N 11	27
70289	 crimp closure N 8	18
702891	 crimp neck vial N 11	27
702891.0	 crimp neck vial N 11	27
702892	 crimp neck vial N 11	26
702893	 screw neck vial N 8	19
	 screw neck vial N 8	34
702894	 pre-sealed vial-closure combination N 24	45
702895	 pre-sealed vial-closure combination N 11	29
702904	 septum N 19	42
702917	 crimp neck vial N 20	38
702918	 crimp neck vial N 20	38
702924	 crimp neck vial N 20	38
702925	 screw closure N 13	33
702926	 screw closure N 13	33
702927	 crimp closure N 20	40
702928	 crimp closure N 20	41
702928.9	 crimp closure N 20	41
702929	 crimp closure N 20	41
702931.1	 stopper N 20	42
702931.3	 stopper N 20	42
702962	 screw neck vial N 13	33
	 screw neck vial N 13	34
702963	 screw cap N 13	33
702966	 screw cap N 13	33
702968	 insert for vials N 8	19
	 insert for vials N 11, small opening	26
702968.1	 insert for vials N 8	19
	 insert for vials N 11, small opening	26
702972	insert for vials N 13	33

Numerical index

REF	Product	Page
702973	 screw neck vial N 13	33
	 screw neck vial N 13	34
702974	 metal spring for insert 702972	33
702981	 septum N 18	37
702995	 crimp closure N 11	28
702995.1	 crimp closure N 11	28
702995.2	 crimp closure N 11	28
702995.3	 crimp closure N 11	28
707971	 linear capillary connector	46
707972	 linear capillary connector	46
707973	 capillary Y connector	46
735111	 standard manual crimper N 11	27
	 standard manual crimper N 11	47
735113	 standard manual crimper N 13	32
	 standard manual crimper N 13	47
735120	 standard manual crimper N 20	39
	 standard manual crimper N 20	47
735126	 standard manual crimper N 8	18
	 standard manual crimper N 8	47
735132	 standard manual crimper N 20	39
	 standard manual crimper N 20	47
735133	 standard manual crimper N 13	32
	 standard manual crimper N 13	47
735208	 ergonomic manual crimper N 8	18
	 ergonomic manual crimper N 8	48
735211	 ergonomic manual crimper N 11	27
	 ergonomic manual crimper N 11	48
735220	 ergonomic manual crimper N 20	39
	 ergonomic manual crimper N 20	48
735311	 ergonomic manual decapper N 11	27
	 ergonomic manual decapper N 11	48
735320	 ergonomic manual decapper N 20	39
	 ergonomic manual decapper N 20	48
735408	 standard manual decapper N 8	18
	 standard manual decapper N 8	47
735500	 battery for electronic crimping tools	48, 49
735501	 stand for electronic crimping tools	49
735509	 rack for crimping tools	48, 49
735511	 electronic crimper N 11	27
	 electronic crimper N 11	48
735513	 electronic crimper N 13	32
	 electronic crimper N 13	48
735520	 electronic crimper N 20	39
	 electronic crimper N 20	48
735532	 electronic crimper N 20	39
	 electronic crimper N 20	48
735533	 electronic crimper N 13	32
	 electronic crimper N 13	48
735611	 electronic decapper N 11	27
	 electronic decapper N 11	48
735613	 electronic decapper N 13	32
	 electronic decapper N 13	48
735620	 electronic decapper N 20	39
	 electronic decapper N 20	48


REF	Product	Page
735700	 high power crimping tool (HP)	27
	 high power crimping tool (HP)	32
	 high power crimping tool (HP)	39
	 high power crimping tool (HP)	49
735708	 crimping head N 8 for HP crimping tool	49
735711	 crimping head N 11 for HP crimping tool	27
	 crimping head N 11 for HP crimping tool	49
735713	 crimping head N 13 for HP crimping tool	32
	 crimping head N 13 for HP crimping tool	49
735720	 crimping head N 20 for HP crimping tool	39
	 crimping head N 20 for HP crimping tool	49
735732	 crimping head N 20 for HP crimping tool	39
	 crimping head N 20 for HP crimping tool	49
735733	 crimping head N 13 for HP crimping tool	32
	 crimping head N 13 for HP crimping tool	49
735811	 decapping head N 11 for HP crimping tool	27
	 decapping head N 11 for HP crimping tool	49
735813	 decapping head N 13 for HP crimping tool	32
	 decapping head N 13 for HP crimping tool	49
735820	 decapping head N 20 for HP crimping tool	39
	 decapping head N 20 for HP crimping tool	49
735911	 standard manual decapper N 11	27
	 standard manual decapper N 11	47
735913	 standard manual decapper N 13	32
	 standard manual decapper N 13	47
735920	 standard manual decapper N 20	39
	 standard manual decapper N 20	47

Tips for correct crimping

For optimal crimp results the appropriate crimping tool needs to be adjustable, as type and height of the vial's crimp neck as well as thickness and hardness of the septum are decisive factors in the crimping process. Also variances from batch to batch within the tolerance of each of the involved components may require adjustments of the tool, in order to achieve a perfect and gas-tight seal.

Basically the following recommendations can be made:

- Do not overcrimp a vial with the intention to make the vial tighter for highly volatile samples. In cases of overcrimped vials the septum is either sucked into the vial or emerging out of the cap's center hole. In both cases the septa material is drawn away from the edges towards the middle of the vial due to the high pressure, leaving the edges insufficiently sealed. Besides that, the convex / concave surface of the septum along with its surface tension may lead to penetration problems (e.g., septum being pushed into the vial by the needle during penetration, coring, etc.).



- Deformation of the cap's sides is also an indication for over-crimping with the same consequences as stated above. Another evidence for over-crimping is a clearly visible formation of wrinkles of the aluminum below the crimp neck.
- In the case of 20 mm pressure release caps the bridges of the score lines may also suffer under too much pressure and thus may burst before reaching the specified level of vial pressure, causing problems in the autosampler.
- Generally beveled top crimp necks N 20 have a lower height than the flat DIN crimp neck. In order to seal beveled top Headspace vials tightly, crimping height and crimping pressure have to be adjusted to a firmer seal.



- Decappers do not need any adjustment.
- To check the crimp result, look for the flatness of the septum and make sure that there are no cap deformations, in order to exclude any signs of over-crimping. To check whether the crimping is too loose, try to turn the cap with 2-3 fingers only. For turning, do not use the whole hand or a large amount of power, because this may lead to an impression of under-crimping which is, however, incorrect. The PTFE lamination of the septum has a very smooth surface, which may slip on the glass rim of the vial allowing to turn the cap already with a relatively low expenditure of energy and the right torque. This is even easier in the case of a beveled top Headspace crimp neck, since on the tapered edge of the rim the septum only has a small area of support.
- Re-usage of crimp neck vials is dangerous, as repeated washing along with the pressure of several crimping and decapping processes may result in neck breakage.



Our program comprises the following types of crimping tools:

- Manual crimpers and decappers (standard) for crimp caps N 8, N 11, N 13 and N 20 as well as for Flip Top/Flip Off caps N 13 and N 20
- Manual crimpers and decappers (ergonomic) for crimp caps N 8 (only crimpers), N 11 and N 20
- Electronic crimpers and decappers (battery-powered) for crimp caps N 11, N 13 and N 20 (not suitable for magnetic and bi-metal caps) as well as for Flip Top / Flip Off caps N 13 and N 20.
- Electronic high power crimping tool with power supply and with exchangeable heads for crimping and decapping of crimp caps N 8 (crimping head only), N 11, N 13 and N 20 as well as crimping heads for Flip Top / Flip Off caps N 13 and N 20.








These products are listed in the chapter of the respective crimp neck diameter as well as in the chapter "Vial accessories and other consumables for chromatography" starting on page 47.

Real size drawings

Screw neck vials N 8 Screw neck vials N 9




Page 19			Page 21									
			Poly-propylene  PP / Glass  Poly-propylene 									
70213	702004	702860	702282	702283	702500	702006	702007	702007.0	702135	702088	702009	702010
70213.2	702893		702293	702284			702008		702335		702172	

Vials N 10 Crimp neck vials N 11

Page 25		Pages 26, 27								
		small opening	wide opening	wide opening						PP / Glass
										
702011	702012	70201CG	70201HP	702885	702888	702141	702015	702891	702891.0	702134
702013		70214CG	70201HP.2	702892			702016	702014		702334


Snap ring vials N 11

Crimp and screw neck vials N 13

Page 30						Page 32		Page 33		
										
702714	702713	702709	702134	702809	702174	70203		702962		702972
702712		702708	702334	702173				702973		





Real size drawings

Crimp neck vials N 8

Page 18		Page 18
		
70286	70282	70202.1





Shell vials N 8

Inserts for vials N 8, N 9, N 10 and N 11

Pages 19, 26 for small opening vials				Pages 21, 25, 26, 30 for wide opening vials				
								
702968	702968.1	702824	702005	702813 702077 s	702716	702818 702818.1 s	702819	702825

Snap cap vials N 18 and N 22


Screw neck vials N 18

Page 36				Page 37		
						
70271	70272	702019	70273	702866 702866.2	702826 702826.2	702826.0

Real size drawings


Crimp neck vials N 20: 5–10 mL

Page 38





Crimp neck vials N 20: 20 mL

Page 38





Real size drawings

Crimp neck vials N 20: 25 mL, 50 mL, 100 mL

Page 39		
		
70210.36	70208.36	70209.1

Screw neck vials for storage of liquid samples

Page 35			
N 15		N 18	N 20
			
702096 / 702311	70285 / 702097	702098	702099

 Screw neck vials N 24



702021 / 702022

702132 / 702133

702023 / 702024

702074 / 702131

General remarks

Translation of cross-references by the VialFinder under www.mn-net.com/vialFinder

The VialFinder is a database-driven translation tool for cross-references of instrument manufacturers and suppliers of consumables worldwide.

The VialFinder immediately shows all options available from MACHEREY-NAGEL for the product of interest. The Finder shows 1:1 matches (in bold type) as well as possible alternative products (in normal type) that – in spite of technical differences to the indicated product – are suitable for the application. The corresponding link on the MN REF will lead you to the appropriate product page on our website that will give information on technical product features as well as possible illustrations of the product.

In case you cannot find your part number via the VialFinder, please send your inquiry by e-mail to support@mn-net.com providing us with all product information you may have. We will then check, if we can offer an equivalent product.



Technical characteristics of vials and closures

Except for the snap cap vials for storage of powdery samples and the blow-molded glass 70209.1, the vials of our program are made from 1st hydrolytic class glass. The dimensions stated in this catalog with respect to vial diameter and height are exact values. Please note that other suppliers often list rounded values (e.g., 12 x 32 mm instead of 11.6 x 32 mm), the actual dimensions are, however, identical due to the required fit in the instru-

ment. Our data concerning the volume are defined realistically usable volumes, not calculated values. For reasons of safety we state rather low values. Here, too, deviations of data of other suppliers may occur, who either use the calculated volume (e.g., 2 mL instead of 1.5 mL) or a defined, realistically usable volume in the upper range (e.g., 1.8 mL instead of 1.5 mL).

Samples for evaluation

Sample packs of all vials and closures can be requested at any time. The sample packs contain 5 pieces of the respective product.



Optimal crimping

For an optimal crimp result the crimping tool needs to be adjusted to:

- Type and height of the vial's crimp neck
- Thickness and hardness of the septa
- Properties of the cap (type, material)

For doing so, please refer to the instruction manual of the individual tool.

Permanent control of the crimp result and thus of the crimping tool settings is necessary.

Incorrect crimping can be recognized by the following features:



Cap deformation



Pulled up edge of the center hole



Strong formation of wrinkles



Convex looking liner



Cap can be turned with only low expenditure of power

General remarks

Website

The new website offers you support in many ways when searching for the optimal product for your application. Use the digital tools described below of our website in your daily laboratory work!

1. Search by product features:

Open the filter bar and select required individual product features from the drop-down menus. Applicable results are displayed to you in tile view.

2. Translation of cross-references:

Follow the link www.mn-net.com/vialfinder.

With the VialFinder part numbers of more than 50 consumables suppliers can be translated.

3. Product information:

Every product on the website has an extensive product detail page with technical data, some even with an additional "Learn more" tab with valuable background information. There is also a download area with relevant product literature.

4. Additional recommendations:

A selection of products with similar product characteristics, such as closures with the same septum, however, with different cap colors, is displayed on the product detail pages under the "Related Products" tab.

In contrast to that, the "Accessories" tab shows products that can be used in conjunction with the product you are looking for, such as vials that can be used in combination with a closure.

5. MN Information center at www.mn-net.com/chroma-news/

Stay up to date on product launches and interesting topics related to chromatography. Please kindly do so by studying the Chroma News, which provide exciting information. Tags make it easier to choose a topic.

6. Supporting videos under

www.mn-net.com/e-training-chromatography

In the e-training section you will find videos that demonstrate practical and theoretical content in an illustrative way. For example, the principle of sample preparation using the so-called SPE is explained there or you can actively follow how to use electronic crimping tools.

7. Featured analytical topics:

Under "Chromatography Service" you will find in the subsection "Featured topics" specific analytical applications, that are reported on in detail. There you also have access to our application database.

8. Conveniently download needed CoCs/CoAs from our website for your own ISO documentation under:

www.mn-net.com/certificate

And finally: purchase our products fast and conveniently in our webshop - your exclusive shopping center.



Learn more about optimal usage of our vials and caps section on our website by watching our video under:
www.mn-net.com/Inserts-Information

Closure selection in GC / HPLC

The choice of the best closure depends on certain features of the instrument (needle type / design, transportation mechanism of the autosampler, etc.) as well as on the requirements of the application (temperature, sensitivity of the analysis, single / multiple injections, etc.) and thus is more complicated and more individual than selection of the correct vial type.

Basically the following recommendations can be made:

- Due to the relatively thick and blunt HPLC needles, only Silicone / PTFE closures, either with or without slit, should be used in combination with them.
- Screw closures N 9 are universally suitable on most autosamplers, convenient in handling and available in a broad selection of different cap colors and septa materials. They fulfill all requirements with regard to tightness and analytical purity for GC as well as for HPLC. Due to the relatively thin septa penetration is safe and easy. Crimp closures N 11 are also universally suitable with regard to autosampler compatibility, however, they are not as safe and convenient in their closing technique as the screw closures N 9.
- Snap ring closures N 11 should only be used in HPLC, as the punctual compacting pressure of the septum against the vial rim by the four pins in the cap does not achieve the same level of tightness as the evenly applied pressure through a circular thread or by crimping.
- For sensitive analyses only high purity Silicone / PTFE closures can be used; if additionally there is a need for minimal coring during penetration, a PTFE / Silicone / PTFE septum (sandwich septum) is recommendable.
- Cap colors may be used for marking (sample marking / lab marking / shift marking). However, please consider that some autosamplers working with photocells may not be able to recognize transparent caps.
- Robotic sample preparation lines as well as some autosamplers are constructed that way that they lift and transport vials with a magnet the size of the vial's cap. In these cases aluminum crimp or plastic screw caps cannot be used, as these materials won't react to a magnet. Only metals like steel or tin will do so. Please check in the catalog for magnetic screw closures N 9, magnetic crimp closures N 11, magnetic screw closures N 18 and magnetic or bi-metal crimp closures N 20.
- For sample storage closed top screw closures (without center hole) should be used. Generally, these also need an elastomeric liner for sealing vials with liquid samples tightly.
- The ability to penetrate a septum depends on the hardness and the total thickness of the septum. The hardness (durometer) is indicated in "° shore". The higher the number, the harder the septum. For easily bending needles as well as for blunt needles with penetration problems soft and preferably thin septa should be chosen. This supports needle safety (clogging/bending of the needle) as well as reliability of the analysis (no contamination of the sample due to particles generated during penetration).
- Due to their artificially reduced cap height screw caps N 9 don't have a standardized thread design. Therefore, it is recommendable only to use vials and closures from one source of supply, in order to ensure a harmonious and tight matching of both components.
- Replacement septa are partially available, however, in case of manual assembly you have the risk of contamination with skin fat / sweat and of a possible wrong side orientation. Therefore we highly recommend only to use ready assembled closures, where the liner perfectly matches the cap and has been automatically inserted under strict hygienic conditions.
- Normally ready assembled closures should be suitable for all types of needles, provided the proper type of septum has been selected. Nevertheless, there might be cases where usage of bonded closures (cap and liner form an inseparable unit) can be recommendable. Example: blunt HPLC needle, however, due to the risk of sample loss / concentration changes no septa with slit can be used. In order to avoid that the unslit septum is pushed into the vial by the needle, you use a bonded closure with unslit septum.
- Almost all liners are covered with PTFE as the chemically resistant barrier. However, when you are doing PFAS analysis, it is crucial to select the right vials and closures for this application. Adsorption effects of glass as well as possible contaminations of the sample by particles from the septa, especially from the PTFE lamination, may put your analysis results at risk. MACHEREY-NAGEL offers polypropylene screw neck vials N 9 and snap ring vials N 11 as well as appropriate closures with a silicone septum laminated either with polyimide or polypropylene. In contrast to PTFE covered liners, these are fluorine-free.



General remarks

Septa guide

Septa material	Temperature resistance from / to	Analytical purity	Fragmentation due to hardness and molecular structure (coring)	Hardness (needle penetration)	Resealability (in case of multiple injections)
PTFE virginal	-200 °C / +260 °C	very high		very hard (but very thin material)	no resealability
Natural rubber / PTFE	-40 °C / +120 °C	low	high, big particles	very hard	high
Red Rubber / TEF (FEP)	-40 °C / +110 °C	medium	medium	medium hard	medium
Butyl	-40 °C / +120 °C	medium	medium	medium hard	medium
Butyl / PTFE	-40 °C / +120 °C	medium	medium	medium hard	medium
Silicone / PTFE	-60 °C / +200 °C	high	low to medium	soft	low to medium
PTFE / Silicone / PTFE	-60 °C / +200 °C	high	very low	soft	very low

Trademarks (registered and common law) of other companies

Waters®: Waters Corp (USA); TurboMatrix™: PerkinElmer® (USA); Merck®: Merck (Germany); CenterGuide™: CRS, USA; Dionex®: Dionex Corporation; Hitachi®: Hitachi, Ltd.; Agilent®: Agilent Technologies, Inc; Shimadzu®: Shimadzu Corporation; Thermo Scientific®: Thermo Fisher Scientific Inc; Varian®: Varian, Inc.; Metrohm: Metrohm AG, Switzerland

Disclaimer

All used names and denotations can be brands, trademarks or registered labels of their respective owner – also if they do not have a special denotation. To mention products and brands is only a kind of information, i.e. it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment.

Regarding these products or services we can not grant any guarantees regarding selection, efficiency or operation.

Product use restriction

MACHEREY-NAGEL chromatography products are intended, developed, designed and sold FOR RESEARCH AND DEVELOPMENT PURPOSES AND ANALYTICAL QUALITY CONTROL / ROUTINE MEASUREMENTS ONLY, except, however, any other function of the product being EXPRESSLY set forth in original MACHEREY-NAGEL product leaflets.

MACHEREY-NAGEL products are intended for GENERAL LABORATORY USE ONLY!

MACHEREY-NAGEL products are suited for QUALIFIED PERSONNEL ONLY !

MACHEREY-NAGEL products shall in any event be used wearing adequate PROTECTIVE CLOTHING!

MACHEREY-NAGEL products shall exclusively be used in an ADEQUATE TEST ENVIRONMENT.

MACHEREY-NAGEL does not assume any responsibility for damages due to improper application of our products in other fields of application.

The user has to ensure that the products used are suitable for the intended application.

General remarks

All information is subject to technical changes. All product data are subject to the currently valid specifications.



Other useful MACHEREY-NAGEL chromatography products from sample preparation to subsequent analysis

SPE



- Concentration and clean-up of your target compounds
- CHROMABOND® and CHROMAFIX® columns and cartridges
- Flash columns, packing material and accessories

Sample filtration



- CHROMAFIL® and CHROMAFIL® Xtra syringe filters
- Protection of sensitive and expensive laboratory equipment
- Use our online FilterFinder on www.mn-net.com/filterfinder

HPLC



- NUCLEOSIL® – the original, one of the first spherical HPLC silicas
- NUCLEODUR® – professional solutions by high purity spherical silica
- NUCLEOSHELL® – highest efficiency with core-shell HPLC columns

TLC



- Quick, easy and cost-efficient separation
- Screening applications and routine analysis
- TLC glass plates, aluminium and polyester sheets as well as accessories

GC




- Highly inert and low-bleed OPTIMA® GC capillary columns
- GC columns for routine analysis and special separations
- More than 50 stationary phases



MACHEREY-NAGEL – Everything you need for successful work in your chromatography lab.

Please also visit our application database with more than 3,000 applications under www.mn-net.com/apps



General remarks

Contacts

With vials and caps of MACHEREY-NAGEL you can trust on safe consumables for reliable analysis results. In combination with numerous other chromatography products from us you have the unique chance to buy out of one source your whole analysis requirements in a cost-effective and time saving way.

Our sales team and our partners worldwide will be more than happy to serve your chromatography needs in the best way. Our expertise and knowledge is at your disposal at any time. You can contact us or our subsidiaries as shown below:

Germany

MACHEREY-NAGEL GmbH & Co. KG
Valencienner Str. 11 · 52355 Düren · Deutschland
Tel.: +49 24 21 969-0
Fax: +49 24 21 969-199
info@mn-net.com

Switzerland

MACHEREY-NAGEL AG
Hirsackerstr. 7 · 4702 Oensingen · Switzerland
Tel.: +41 62 388 55 00
Fax: +41 62 388 55 05
sales-ch@mn-net.com

France

MACHEREY-NAGEL SARL à associé unique
1, rue Gutenberg · 67722 Hoerdt · France
Tél. : +33 388 68 22 68
Fax : +33 388 51 76 88
sales-fr@mn-net.com

USA

MACHEREY-NAGEL Inc.
924 Macron Blvd., Suite 102 · Allentown, PA 18109 · USA
Tel.: +1 484 821 0984
Tel.: +1 888 321 6224 toll free
sales-us@mn-net.com

An overview of dealers who sell our products worldwide can be found on our website at www.mn-net.com/distributors.

For specific technical questions on vials and caps, please contact:

Phone: +49 24 21 969-333
e-mail: support@mn-net.com

Your local distributor

KATEN200010 ValsCat en6 / 12.2025

www.mn-net.com

MACHEREY-NAGEL



Management
System
EN ISO 13485:2016
ISO 9001:2015

www.tuv.com
ID 0000056401



MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com
US +1 888 321 62 24 sales-us@mn-net.com