

**Data sheet**

08.02.2023

# Semi-quantitative test strips QUANTOFIX Potassium

Parameter*	Value*
REF	91316
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Potassium
Gradation	0 • 200 • 400 • 700 • 1000 • 1500 mg/L K <sup>+</sup>
Method	Dipikrylamine
Color change	yellow → orange

\*These data are typical values and must not be regarded as specification.

