

## Data sheet

08.02.2023

## Semi-quantitative test strips QUANTOFIX Nitrite

Parameter*	Value*
REF	91311
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Nitrite
Gradation	0 • 1 • 5 • 10 • 20 • 40 • 80 mg/L NO <sub>2</sub> <sup>-</sup>
Method	Griess reaction
Color change	white → red-violet
Instrumental measuring range	0.5–80 mg/L NO <sub>2</sub> <sup>-</sup>

\*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11  
52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com