

Data sheet

Filter paper

Parameter*	Value*
Grade	MN GF-1
Basis weight	52 g/m ²
Thickness	0.28 mm / 11 mil
Filtration speed	12 s
Maximum temperature	500 °C / 932 °F
Surface	Smooth
Color	White
Material	Borosilicate glass micro fibers
Average retention capacity	0.7 µm
Characteristics / Application	Glass fiber filter without binder, fast filtration and very good retention capacity, for high efficiency laboratory filtration and air sampling

*These data are typical values and must not be regarded as specification.