






visocolor® ECO Zink
 Zinc / Zinc / Zinc

 Messbereich / Range / Domaine de mesure / Rango
 0.5 – 3 mg/L Zn²⁺

 Probe
 Sample
 Echantillon
 Muestra





2 ∆ Zn-1


 Schütteln
 Shake
 Agiter
 Agitar


2 ∆ Zn-2


 Schütteln
 Shake
 Agiter
 Agitar


5 ∆ Zn-3


 Schütteln
 Shake
 Agiter
 Agitar

 Warten
 Wait
 Attendre
 Espere

 Messung
 Measurement
 Mesure
 Medición

 Geeignet für Meerwasser nach Verdünnung 1 + 9
 Suitable for sea water after dilution 1 + 9
 Convient à l'eau de mer après dilution 1 + 9
 Adecuado para el agua de mar tras su dilución 1 + 9