

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 85/90 BF
Basis weight	90 g/m ²
Thickness	0.51 mm / 20.1 mil
Filtration speed	15 s
Maximum temperature	500 °C / 932 °F
Surface	Smooth
Color	White
Material	Borosilicate glass micro fibers
Average retention capacity	0.5 µm
Characteristics / Application	Glass fiber filter without binder, suitable for determination of suspended solids according to DIN EN 872:2005

*These data are typical values and must not be regarded as specification.