

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Total iron 100

Parameter*	Value*
REF	91344
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Iron
Gradation	0 • 2 • 5 • 10 • 25 • 50 • 100 mg/L Fe ²⁺ 3 ⁺
Method	Triazine
Color change	white → blue-violet

*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG
Valenciener Str. 11
52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
US Tel.: +1 888 321 62 24 sales-us@mn-net.com