

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Nitrate / Nitrite

| Parameter* | Value* |
|------------------------------|--|
| REF | 91313 |
| Evaluable on reflectometer | QUANTOFIX Relax |
| Shelf life (from production) | 2.5 Years |
| Scope of delivery | 100 test strips in a tube, instruction leaflet |
| Parameter | Nitrate |
| Gradation | 0 • 10 • 25 • 50 • 100 • 250 • 500 mg/L NO ₃ ⁻ |
| Method | Modified Griess reaction |
| Color change | white → red-violet |
| Instrumental measuring range | 10–500 mg/L NO ₃ ⁻ |
| Parameter 2 | Nitrite |
| Gradation | 0 • 10 • 25 • 50 • 100 • 250 • 500 mg/L NO ₂ ⁻ |
| Method | Griess reaction |
| Color change | white → red-violet |
| Instrumental measuring range | 0.5–50 mg/L NO ₂ ⁻ |

*These data are typical values and must not be regarded as specification.

