

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 616 LSA-50
Maximum temperature	100 °C / 212 °F
Surface	Smooth
Color	White
Material	Cellulose, cation exchanger
Capacity	0,05 ± 0,01 meq/cm ² , 0,9 ± 0,2 meq/Filter with ø 48 mm
Characteristics / Application	Filter paper with strongly acidic cation exchange resin, Matrix polystyrene crosslinked with DVB, functional groups SO ₃ H, supplied in H ⁺ form/

*These data are typical values and must not be regarded as specification.