

## Data sheet

### Filter paper

| Parameter*                    | Value*                                                                             |
|-------------------------------|------------------------------------------------------------------------------------|
| Grade                         | MN GF-4                                                                            |
| Basis weight                  | 120 g/m <sup>2</sup>                                                               |
| Thickness                     | 0.53 mm / 20.9 mil                                                                 |
| Filtration speed              | 5 s                                                                                |
| Maximum temperature           | 500 °C / 932 °F                                                                    |
| Surface                       | Smooth                                                                             |
| Color                         | White                                                                              |
| Material                      | Borosilicate glass micro fibers                                                    |
| Average retention capacity    | 1.4 µm                                                                             |
| Characteristics / Application | Glass fiber filter without binder, used e.g. as effective pre-filter for membranes |

\*These data are typical values and must not be regarded as specification.