

Data sheet

Filter paper

| Parameter* | Value* |
|-------------------------------|--|
| Grade | MN 1672 |
| Basis weight | 85 g/m ² |
| Thickness | 0.13 mm / 5.1 mil |
| Filtration speed | 35 s |
| Surface | Smooth |
| Color | White |
| Material | Cellulose |
| Ash content | ca. 0.1 % |
| Average retention capacity | 4–7 µm |
| Characteristics / Application | Wet-strengthened, medium fast filtration |

*These data are typical values and must not be regarded as specification.