

## Data sheet

### Filter paper

| Parameter*                    | Value*                                             |
|-------------------------------|----------------------------------------------------|
| Grade                         | MN 522                                             |
| Basis weight                  | 80 g/m <sup>2</sup>                                |
| Thickness                     | 0.34 mm / 13.4 mil                                 |
| Filtration speed              | 22 s                                               |
| Surface                       | Creped                                             |
| Color                         | White                                              |
| Material                      | Cellulose                                          |
| Ash content                   | ca. 0.1 %                                          |
| Characteristics / Application | Medium fast filtration, for general laboratory use |

\*These data are typical values and must not be regarded as specification.