

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Copper

Parameter*	Value*
REF	91304
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Copper
Gradation	0 • 10 • 30 • 100 • 300 mg/L Cu ²⁺
Method	Biquinoline
Color change	white → red-violet

*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG
 Valencienner Str. 11
 52355 Düren · Germany
www.mn-net.com

DE Tel.: +49 24 21 969-0 info@mn-net.com
 CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com
 FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com
 US Tel.: +1 888 321 62 24 sales-us@mn-net.com