

Customer information

Version 5

List of functions accessible and blocked in user mode on the VIS II and UV-VIS II

The following lists contains the functions that are active or blocked in user mode.

Functions	Active	Blocked
Menu 'Dataexport'		
ACRON settings		x
LIMS settings		x
Network Export settings		x
Menu 'Settings'		
Language	✓	
Region	✓	
Date & time		x
Memory settings		x
Signal loudness	✓	
Printer settings	✓	
Remember information		x
Reaction time		x
Lamp control		x
NTU check		x
Sipper settings		x
Backup		x
Dilution format	✓	
Energy mode	✓	
Menu 'System'		
Update		x
Reset		x
Menu 'Calibration'		
Zero calibration		x
Turbidity calibration		x
Remove turbidity values and save changes		x
Add turbidity values and save changes		x

Menu 'User accounts'		
Edit own account	✓	
Update own account to admin		x
Edit another user account		x
Delete a user		x
Create a new user		x
Terminate user mode		x
Measurements		
Delete own measurements		x
Delete other measurements		x
User defined methods		
Create new method	✓	
Export own method	✓	
Export other methods	✓	
Delete methods		x
Edit methods		x

Contact

If you have further questions, please do not hesitate to contact us:

Technical support and customer service

Tel.: +49 24 21 969-333
E-Mail: support@mn-net.com

For orders (international)

E-Mail: sales@mn-net.com