



# visicolor® ECO Nitrit

Nitrite / Nitrite / Nitrito

Messbereich / Range / Domaine de mesure / Rango  
0.02 – 0.5 mg/L NO<sub>2</sub><sup>-</sup>




1  
Probe  
Sample  
Echantillon  
Muestra




2  
4 δ NO<sub>2</sub><sup>-</sup>1




3  
Schütteln  
Shake  
Agiter  
Agitar




4  
1 → NO<sub>2</sub><sup>-</sup>2



5  
Schütteln  
Shake  
Agiter  
Agitar



6  
Warten  
Wait  
Attendre  
Espere



7  
Messung  
Measurement  
Mesure  
Medición

Geeignet für Meerwasser  
Suitable for seawater  
Convient à l'eau de mer  
Adecuado para el agua de mar

