

NucleoMag® NGS Clean-up (DNA/RNA) – MagnetaPure 32+

Protocol details

Application	PCR Clean-up
Kit	NucleoMag® NGS Clean-up and Size Select (Protocol 5.4)
REF	744970 (.5 / .50 / .500)
Protocol name	NMNGSCleanup
Run time	approx. 31 min

Rev 01



Four easy steps

Procedure	
1	Fill the 96-well Deep-well plates according to the table below.
2	Load the plates on the MagnetaPure 32+.
3	Insert tip combs on the mounting grooves.
4	Select the protocol from the instrument menu and start run.

Note: Please equip all tip combs in order to cover the magnetic rods in used and unused wells.

Loading table for plate 1 and 2

Position	Reagents	Samples per plate
Column 1 + 7	DNA / RNA sample* (100 µL), NucleoMag® NGS Beads (180 µL)	Sample 1-8 Sample 9-16
Column 2 + 8	80% Ethanol (200 µL)	Sample 1-8 Sample 9-16
Column 3 + 9	80% Ethanol (200 µL)	Sample 1-8 Sample 9-16
Column 4 + 10	empty	- -
Column 5 + 11	empty	- -
Column 6 + 12	Elution Buffer (25 µL)**	Sample 1-8 Sample 9-16

Please refer to the image below for a visual representation of the loading scheme

Note: * Sample volume can be adjusted (10 -100 µL); Ratio between sample volume and NGS Bead Suspensions should be kept constant (1:1.8)

** Elution can be performed in a minimal volume of 20 µL

Loading scheme

	1	2	3	4	5	6	7	8	9	10	11	12	
A													
B													
C													
D													
E													
F													
G													
H													

Additional consumables and instrumentation

Product	Specification	REF
MagnetaPure 32+	Automated nucleic extraction system for MACHEREY-NAGEL's NucleoMag® kits enabling parallel processing of up to 32 samples	747010
96 Deep-well Plate	96 deep-well plates for MagnetaPure 32+ (25 pieces)	744955
Tip Combs	8-well tip combs for MagnetaPure 32+ (50 pieces)	744960

Disclaimer

Information

MACHEREY-NAGEL GmbH & Co. KG makes every effort to include accurate and up-to-date information within this publication; however, it is possible that omissions or errors might have occurred. MACHEREY-NAGEL GmbH & Co. KG cannot, therefore, make any representations or warranties, expressed or implied, as to the accuracy or completeness of the information provided in this publication. Changes in this publication can be made at any time without notice. For technical details and detailed procedures of the specifications provided in this document please contact your MACHEREY-NAGEL representative. This publication may contain reference to applications and products which are not available in all markets. Please check with your local sales representative. All mentioned trademarks are protected by law. All used names and denotations can be brands, trademarks, or registered labels of their respective owner – also if they are not special denotation. To mention products and brands is only a kind of information (i.e., it does not offend against trademarks and brands and can not be seen as a kind of recommendation or assessment). Reference herein to any specific commercial product, process, or service by trade name, trademark, manufacturer, or otherwise, does not necessarily constitute or imply its endorsement, recommendation, or support by MACHEREY-NAGEL GmbH & Co. KG. Any views or opinions expressed herein by the authors' do not necessarily state or reflect those of MACHEREY-NAGEL. NucleoMag® is a registered trademark of MACHEREY-NAGEL GmbH & Co. KG, Düren, Germany.