


NANOCOLOR[®] gesamt-Stickstoff TN_b 60

total Nitrogen / Azote total / Nitrógeno total










Test 0-92
REF 985092

Methode / Method / Méthode / Método
0921 3–60 mg/L N


345 / 350 / 365 nm

Min	20–80 %	Max
3	14 49	60 mg/L N



NANOCOLOR[®] gesamt-Stickstoff TN_b 60

total Nitrogen / Azote total / Nitrógeno total



Test 0-92
REF 985092

Methode / Method / Méthode / Método
0921 3–60 mg/L N

345 / 350 / 365 nm

Min	20–80 %		Max
3	14	49	60 mg/L N

