

## Data sheet

## Filter paper

| Parameter*                    | Value*                                                 |
|-------------------------------|--------------------------------------------------------|
| Grade                         | MN 321                                                 |
| Basis weight                  | 85 g/m <sup>2</sup>                                    |
| Thickness                     | 0.23 mm / 9.1 mil                                      |
| Filtration speed              | 5 s                                                    |
| Surface                       | Smooth                                                 |
| Color                         | White                                                  |
| Material                      | Cellulose                                              |
| Ash content                   | ca. 0.1 %                                              |
| Average retention capacity    | 7–12 µm                                                |
| Characteristics / Application | Nitrogen-free, fast filtration, e.g. for beer analysis |

\*These data are typical values and must not be regarded as specification.