

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 640 md
Basis weight	85 g/m ²
Thickness	0.18 mm / 7.1 mil
Filtration speed	55 s
Surface	Smooth
Color	White
Material	Cellulose
Ash content	< 0.01 %
Average retention capacity	4–7 µm
Characteristics / Application	Ashless, medium fast filtration, e.g. for magnesium ammonium phosphate, zinc sulphide, etc.

*These data are typical values and must not be regarded as specification.