

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Chlorine

Parameter*	Value*
REF	91317
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Chlorine
Gradation	0 • 1 • 3 • 10 • 30 • 100 mg/L Cl ₂
Method	Redox reaction
Color change	white → red-violet

*These data are typical values and must not be regarded as specification.

