

**Data sheet**

08.02.2023

# Semi-quantitative test strips QUANTOFIX Total sugar

Parameter*	Value*
REF	91352
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Fructose / Glucose
Gradation	0 • 55 • 100 • 250 • 400 • 600 • 800 mg/L Fructose / Glucose
Method	Enzymatic
Color change	yellow → ocher
Instrumental measuring range	55–700 mg/L Fructose / Glucose

\*These data are typical values and must not be regarded as specification.

