



visicolor® ECO Chrom(VI)

Chromium(VI) / Chrome(VI) / Cromo(VI)

Messbereich / Range / Domaine de mesure / Rango
0.02 – 0.50 mg/L Cr(VI)




1
Probe
Sample
Echantillon
Muestra




2
5 ⚡ Cr-1



3
Schütteln
Shake
Agiter
Agitar




4
5 ⚡ Cr-2



5
Schütteln
Shake
Agiter
Agitar



6
Warten
Wait
Attendre
Espere



7
Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar

