

## Data sheet

08.02.2023

## Semi-quantitative test strips QUANTOFIX Chlorine dioxide 15

Parameter*	Value*
REF	91355
Evaluable on reflectometer	No
Shelf life (from production)	2 Years
Scope of delivery	50 test strips in a tube, instruction leaflet
Parameter	Chlorine dioxide
Gradation	0.25 • 0.5 • 1 • 3 • 5 • 10 • 15 mg/L ClO <sub>2</sub>
Method	Redox reaction
Color change	colorless → pink

\*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11  
52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com