

Data sheet

Filter paper

Parameter*	Value*
Grade	MN 85/70 BF
Basis weight	70 g/m ²
Thickness	0.33 mm / 13.0 mil
Filtration speed	15 s
Migration distance	135–165 mm/10 min
Maximum temperature	500 °C / 932 °F
Surface	Smooth
Color	White
Material	Borosilicate glass micro fibers
Average retention capacity	0.6 µm
Characteristics / Application	Glass fiber filter without binder, free of cellulose

*These data are typical values and must not be regarded as specification.