

MACHEREY-NAGEL

Papers for seed germination testing

Filtration



Premium quality, reliability and flexibility

- Comply with ISTA specifications
- Tailor-made customer solutions
- Made in Germany

MACHEREY-NAGEL

www.mn-net.com



Papers for seed germination tests

ISTA-compliant papers for seed testing

The grades specially selected for seed testing are made of high-purity cellulose and are optimally adapted by MACHEREY-NAGEL's experts to meet our customers' requirements. International seed producers trust in the quality of MN for years.

The papers are characterized by the following properties:

- Made from high-purity cellulose
- Excellent water absorption
- Roots do not grow into the paper
- Conductivity < 40 mS/m
- pH 6.0 – 7.5

MN offers you a wide selection of papers for the ISTA methods:

- Top of paper (TP)
- Pleated paper (PP)
- Between Paper (BP)
- Heubach Testing, ESA-recommended

Your MACHEREY-NAGEL support specialists will guide you from the sampling process, through day-to-day business, to advice on new challenges.

Ordering information

Method	Grade	Description	REF
Pleated paper			
	MN 7229	Pleated strips, 110 × 20 × 2000 mm 50 double pleats, 1008 strips/pack	1007229001
	MN 7232		1007232001
	MN 614	Wrapping strips, embossed	On request *
Top of paper			
	MN 7228	Smooth filter paper, 110 g/m ² , 0.18 mm thickness	On request *
	MN 270	Smooth filter paper, 270 g/m ² , 0.54 mm thickness	On request *
Between paper			
	MN 7230	Creped filter paper, 160 g/m ² , 0.69 mm thickness	On request *
	MN 651/120	Creped filter paper, 120 g/m ² , 0.44 mm thickness	On request *
	MN 7231	Smooth filter paper, 150 g/m ² , 0.32 mm thickness	On request *
Heubach Testing – Measurement of dust reference values			
	MN 85/70 BF	Glass fiber filters, circular, with binder, 50 mm	4040050
Determination of fats			
	MN 645	Extraction thimbles for Soxhlet – also quick extractors	On request *

* Custom dimensions available on request. Please contact us!

If you have any questions or need a quote, please contact us:
E-mail: support@mn-net.com
Phone: 024 21 969-333

www.mn-net.com

MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com

CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com

US +1 888 321 62 24 sales-us@mn-net.com

Follow us



sun_apple (S. 2) – stock.adobe.com
KATENT00262 Flyer Saatgutpapiere en1a

