


NANOCOLOR® gesamt-Phosphat 45

total Phosphate/Phosphate total/Fosfato total













Test 0-55
REF 985055

Methode / Method / Méthode / Método

690 / 710 nm


0551 5.0–50.0 mg/L P
0552 15–150 mg/L PO₄³⁻
0553 10–115 mg/L P₂O₅

Min	20–80 %		Max
5.0	14.0	41.0	50.0 mg/L P
15	42	123	150 mg/L PO ₄ ³⁻
10	31	94	115 mg/L P ₂ O ₅



NANOCOLOR[®] ortho-Phosphat 45

ortho-Phosphate / ortho Phosphates / orto Fosfato



Test 0-55
REF 985055

Methode / Method / Méthode / Método

690 / 710 nm

0554 5.0 – 50.0 mg/L PO₄-P
0555 15 – 150 mg/L PO₄³⁻
0556 10 – 115 mg/L P₂O₅

Min	20 – 80 %		Max
5.0	14.0	41.0	50.0 mg/L P
15	42	123	150 mg/L PO ₄ ³⁻
10	31	94	115 mg/L P ₂ O ₅

