

## Data sheet

## Filter paper

Parameter*	Value*
Grade	MN 126/70
Basis weight	70 g/m <sup>2</sup>
Thickness	0.2 mm / 7.9 mil
Filtration speed	25 s
Surface	Creped
Color	White
Material	Cellulose
Characteristics / Application	Medium fast filtration, wet strengthened, for technical applications. Both the pulps used and the wet-strengthening agent are approved under German foodstuffs legislation (Recommendation XXXVI of the BgVV).

\*These data are typical values and must not be regarded as specification.