

Data sheet

Filter paper

Parameter*	Value*
Grade	MN GF-4
Basis weight	120 g/m ²
Thickness	0.53 mm / 20.9 mil
Filtration speed	5 s
Maximum temperature	500 °C / 932 °F
Surface	Smooth
Color	White
Material	Borosilicate glass micro fibers
Average retention capacity	1.4 µm
Characteristics / Application	Glass fiber filter without binder, used e.g. as effective pre-filter for membranes

*These data are typical values and must not be regarded as specification.