



visocolor® ECO Nitrit

Nitrite / Nitrite / Nitrito

Messbereich / Range / Domaine de mesure / Rango

0.02 – 0.5 mg/L NO_2^- 


Probe
Sample
Echantillon
Muestra




4 NO_2^- -1




Schütteln
Shake
Agiter
Agitar




1 NO_2^- -2



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser

Suitable for seawater

Convient à l'eau de mer

Adecuado para el agua de mar