


NANOCOLOR® Ammonium 2000

Amonio



Test 0-02
REF 985002

Methode / Method / Méthode / Método

0021 300 – 1600 mg/L $\text{NH}_4\text{-N}$

0022 400 – 2000 mg/L NH_4^+

0023 400 – 2000 mg/L NH_3

585 nm

Min	20 – 80 %		Max
300	560	1340	1600 mg/L $\text{NH}_4\text{-N}$
400	720	1680	2000 mg/L $\text{NH}_4^+ / \text{NH}_3$

