


# NANOCOLOR® Organische Säuren 3000

Organic Acids 3000 / Acides organiques 3000 / Ácidos orgánicos 3000



Test 0-50

REF 985050









Methode / Method / Méthode / Método

470 / 490 nm

0501 30–3000 mg/L


0502 0.5–50 mmol/L

| Min  | 20–80 % |      | Max       |
|------|---------|------|-----------|
| 30   | 624     | 2406 | 3000 mg/L |
| 0.50 | 10.4    | 40.1 | 50 mmol/L |



# NANOCOLOR® Organische Säuren 3000

Organic Acids 3000 / Acides organiques 3000 / Ácidos orgánicos 3000



Test 0-50

REF 985050

Methode / Method / Méthode / Método

470 / 490 nm

0501 30–3000 mg/L

0502 0.5–50 mmol/L

Min

20–80 %

Max

30

624

2406

3000 mg/L

0.50

10.4

40.1

50 mmol/L

