


NANOCOLOR® gesamt-Stickstoff TN_b 22

total Nitrogen / Azote total / Nitrógeno total









Test 0-83
REF 985083

Methode / Method / Méthode / Método
0831 0.5–22.0 mg/L N


345 / 350 / 365 nm

Min	20–80 %		Max
0.5	4.8	17.7	22.0 mg/L N



NANOCOLOR® gesamt-Stickstoff TN_b 22

total Nitrogen / Azote total / Nitrógeno total



Test 0-83
REF 985083

Methode / Method / Méthode / Método
0831 0.5–22.0 mg/L N

345 / 350 / 365 nm

Min	20–80 %		Max
0.5	4.8	17.7	22.0 mg/L N

