



visocolor® ECO Kupfer

Copper / Cuivre / Cobre

Messbereich / Range / Domaine de mesure / Rango
0.1 – 1.5 mg/L Cu²⁺




Probe
Sample
Echantillon
Muestra




5 ∆ Cu-1




Schütteln
Shake
Agiter
Agitar




5 ∆ Cu-2



Schütteln
Shake
Agiter
Agitar



Warten
Wait
Attendre
Espere



Messung
Measurement
Mesure
Medición

Geeignet für Meerwasser
Suitable for seawater
Convient à l'eau de mer
Adecuado para el agua de mar