

Data sheet

08.02.2023

# Semi-quantitative test strips QUANTOFIX Ammonium

Parameter*	Value*
REF	91315
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Ammonium
Gradation	0 • 10 • 25 • 50 • 100 • 200 • 400 mg/L NH <sub>4</sub> <sup>+</sup>
Method	Nessler
Color change	bright yellow → orange
Instrumental measuring range	10–350 mg/L NH <sub>4</sub> <sup>+</sup>

\*These data are typical values and must not be regarded as specification.

