

MACHEREY-NAGEL

# Rapid Tests for Science

Rapid Tests



## Dip & Read Tests

- Quick & easy
- High informative value
- Accurate results



**MACHEREY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)










# Dip & Read Tests for Science

Dip & Read tests allow the easy detection of important parameters directly at the point of interest. They provide reliable results within a few seconds and allow a quick decision on further activities. MN test strips are the ideal products for everyday lab use and are successfully used in science labs around the globe.



## Product information

Test	Product	Application	Gradation	REF
	pH-Fix 0–14	<ul style="list-style-type: none"> <li>Always ready for use</li> <li>Patented technology prevents bleeding of the indicator</li> <li>Brilliant colour charts for exact results</li> </ul>	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	92110
	pH-Fix 0–14 PT	<ul style="list-style-type: none"> <li>pH-Fix tests also available in PlopTop tubes</li> <li>Virtually unbreakable pack</li> <li>Easy to open and close</li> </ul>	pH 0 · 1 · 2 · 3 · 4 · 5 · 6 · 7 · 8 · 9 · 10 · 11 · 12 · 13 · 14	92111
	QUANTOFIX® Peroxide 100	<ul style="list-style-type: none"> <li>Peroxide containing solvents may explode when heated</li> <li>Test allows to check for peroxide in a few seconds</li> <li>Easy Dip &amp; Read procedure for best usability</li> </ul>	0 · 1 · 3 · 10 · 30 · 100 mg/L H <sub>2</sub> O <sub>2</sub>	91312
	Potassium iodide starch paper	<ul style="list-style-type: none"> <li>Detects strong oxidisers (i. e. free chlorine)</li> <li>Easy control of diazotization reactions (excess of nitrite)</li> </ul>	Detection limit 1 mg/L NO <sub>2</sub> <sup>-</sup> or free Cl <sub>2</sub>	90754
	Fluoride Test Paper	<ul style="list-style-type: none"> <li>Detects dangerous hydrofluoric acid without additional reagents</li> <li>Easy control of puddles with unknown liquids</li> <li>Provides immediate results</li> </ul>	Detection limit 20 mg/L F <sup>-</sup>	90750
	Watesmo	<ul style="list-style-type: none"> <li>Test paper for detection of water</li> <li>Reaction colour irreversible</li> <li>Easy control whether solvents are anhydrous</li> </ul>	Detection of water droplets	90609



Follow us



© Leah-Anne Thompson, Tom, NIDeLander · stock.adobe.com  
KATEN100100 Science en2

[www.mn-net.com](http://www.mn-net.com)

# MACHEREY-NAGEL



Management System  
EN ISO 13485:2016  
ISO 9001:2015  
www.tuv.com  
ID 900008338

MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11  
52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com

