

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Sulfite

Parameter*	Value*
REF	91306
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Sulfite
Gradation	0 • 10 • 25 • 50 • 100 • 250 • 500 • 1000 mg/L SO ₃ ²⁻
Method	Nitroprussid / Zn-hexacyanoferrate
Color change	white → salmon
Instrumental measuring range	10–500 mg/L SO ₃ ²⁻

*These data are typical values and must not be regarded as specification.

