

## Data sheet

08.02.2023

## Semi-quantitative test strips QUANTOFIX Arsenic 10

Parameter*	Value*
REF	91334
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, additional reagents, instruction leaflet in a box
Parameter	Arsenic
Gradation	0•0.01•0.025•0.05•0.1•0.5 mg/L As <sup>3+</sup> 5 <sup>+</sup>
Method	Modified Gutzeit test
Color change	white → yellow-brown

\*These data are typical values and must not be regarded as specification.



MACHEREY-NAGEL GmbH & Co. KG  
Valenciener Str. 11  
52355 Düren · Germany  
[www.mn-net.com](http://www.mn-net.com)

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com