

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Sulfate

Parameter*	Value*
REF	91329
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Sulfate
Gradation	< 200 • > 400 • > 800 • > 1200 • > 1600 mg/L SO ₄ ²⁻
Method	Ba-thorine-complex
Color change	red → yellow

*These data are typical values and must not be regarded as specification.

