









visocolor® ECO DEHA
 DEHA / DEHA / DEHA

 Messbereich / Range / Domaine de mesure / Rango
 0.01 – 0.30 mg/L DEHA

 filtrieren wenn trüb
 filter if cloudy
 filtrer si trouble
 filtra si turbido

 Probe
 Sample
 Echantillon
 Muestra


4 ∅ DEHA-1


 Schütteln
 Shake
 Agiter
 Agitar


4 ∅ DEHA-2


 Schütteln
 Shake
 Agiter
 Agitar

 Warten
 Wait
 Attendre
 Espere

 Messung
 Measurement
 Mesure
 Medición

 Geeignet für Meerwasser
 Suitable for seawater
 Convient à l'eau de mer
 Adecuado para el agua de mar