

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Glucose

Parameter*	Value*
REF	91348
Evaluable on reflectometer	QUANTOFIX Relax
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Glucose
Gradation	0 • 50 • 100 • 250 • 500 • 1000 • 2000 mg/L Glucose
Method	Enzymatic
Color change	yellow → blue-green
Instrumental measuring range	50–2000 mg/L Glucose

*These data are typical values and must not be regarded as specification.

