

Data sheet

08.02.2023

Semi-quantitative test strips QUANTOFIX Nitrite / pH

Parameter*	Value*
REF	91338
Evaluable on reflectometer	No
Shelf life (from production)	2.5 Years
Scope of delivery	100 test strips in a tube, instruction leaflet
Parameter	Nitrite
Gradation	0 • 1 • 5 • 10 • 20 • 40 • 80 mg/L NO ₂ ⁻
Method	Griess reaction
Color change	white → red-violet
Parameter 2	pH
Gradation	6.0 • 6.4 • 6.7 • 7. • 7.3 • 7.6 • 7.9 • 8.2 • 8.4 • 8.6 • 8.8 • 9.0 • 9.3 • 9.6
Method	mixed indicator
Color change	yellow-orange → violet-red

*These data are typical values and must not be regarded as specification.

